

**Southeast  
Medical  
Products, Inc**

*Specializing in  
Orthopedics*

*Since 1977*

*2008 Product Catalog*

For over 30 years, Southeast Medical Products has provided the professional medical community with a great selection of quality products at fair and competitive prices.

We began in 1977 with a single product line, and today we represent numerous manufacturers and warehouse thousands of items. This allows us to offer a wide range of products, and have the ability to ship most orders the same day.

Southeast Medical Products prides itself on providing great service.

Every customer has different needs and we always strive to do whatever it takes to take care of them.

*WE SINCERELY THANK YOU FOR YOUR BUSINESS!*

# OUR MANUFACTURERS...



GALVESTON MEDICAL MFG.



...and many more.

Southeast Medical Products, Inc

5524 Commerce Dr

Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)

800.330.0890 (ph) • 800.851.9012 (fax)

# ACROMIOCLAVICULAR BRACE

Sturdy cotton/duck construction and felt shoulder straps provide additional comfort on shoulder area. Splint helps stabilize shoulder, humerus, and elbow following injuries. One size fits most patients.



Part Number:	Circ.	Size
79-85010	21"-42"	U

**Measurement Information:**

Chest circumference measurement taken at bottom of sternum.


Southeast Medical Products, Inc  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# BOXER'S FRACTURE SPLINTS

## TKO™ (The Knuckle Orthosis)



- Indicated for ligament injuries, Boxer's Fracture and other fractures of the metacarpals and proximal phalanges
- Straps are infinitely adjustable for fit and comfort (pain & swelling)
- Distal Buddy Strapping eliminates rotation of the finger
- Padded stays are malleable
- Reversible - can be used as a radial gutter splint



## ComfortFORM™ Boxer's Splint



Durable lightweight fabric is **Lycra®** lined for breathability and patient comfort. The **ComfortFORM** Boxer's Splint is designed to eliminate the need for ulnar gutter casting or splinting. The malleable ulnar stay creates a dorsal block that maintains the 4th & 5th phalange at a 70° to 90° angle to keep fingers in full flexion following injury. The preformed palmer stay helps provide correct fit and proper support. Indications include Boxer's Fractures, flexor tendon and/or digital nerve damage to the 4th & 5th metacarpal phalange.



## Boxer Splint

Anatomically designed to provide comfortable support, protection and immobilization of the 4th and 5th metacarpals

Indicated for treatment of acute or sub-acute fractures of the 4th and 5th metacarpals commonly referred to as "boxer fractures"

Constructed of lightweight low-density polyethylene shell with soft foam liner and palmer pad and contact closure straps



**Extended Length**

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# CANDO® EXERCISE BANDS

*Exercise band therapy is the most commonly prescribed exercise modality for rehabilitation, conditioning and training. Simply cut an appropriate length of band and begin your exercise. The exercise band can be used with any band or tubing accessory.*



Low Powder Latex Bands		
	<u>50 yd roll, 1/box</u>	<u>6 yd roll, 1/box</u>
	<u>order #</u>	<u>order #</u>
Tan	10-5220	10-5210
Yellow	10-5221	10-5211
Red	10-5222	10-5212
Green	10-5223	10-5213
Blue	10-5224	10-5214
Black	10-5225	10-5215
Silver	10-5226	10-5216
Gold	10-5227	10-5217
Latex Free Bands		
	<u>50 yd roll, 1/box</u>	<u>6 yd roll, 1/box</u>
	<u>order #</u>	<u>order #</u>
Tan	10-5620	10-5610
Yellow	10-5621	10-5611
Red	10-5622	10-5612
Green	10-5623	10-5613
Blue	10-5624	10-5614
Black	10-5625	10-5615
Silver	10-5626	10-5616
Gold	10-5627	10-5617

Clinic and Hospital Packs - Bulk rolls are packaged in 6 yd (5.5 m) and 50 yd (45.7 m) dispenser boxes. The boxes are color-coded to resistance level and design-coded to indicate whether the box contains low-powder or no-latex Cando exercise band, also helpful for inventory control and storage.



**Southeast Medical Products, Inc**  
**5524 Commerce Dr**  
**Orlando FL 32839**

**407.851.5848 (ph) • 407.851.9012 (fax)**  
**800.330.0890 (ph) • 800.851.9012 (fax)**

# CASTING MATERIALS

Whatever brand you prefer, you can count on ***Southeast Medical Products*** to deliver great pricing, quick delivery & unparalleled service.



## 3M™ Scotchcast™ Plus Casting Tape

- Water-activated resin lets you control working time according to application needs.
- High and early strength allows for weight bearing in 20 minutes and provides a durable cast with less chance of breakdown and the need for costly repairs or replacements.
- Lightweight for easy mobility.
- Stretches to conform well to contours for a more comfortable, better fitting cast.
- Porous weave allows the skin to breathe.
- Available in 11 colors to maximize patient satisfaction.



## ***Delta-Lite® Conformable Fiberglass Cast Tape***



- ✓ Exceptional conformability
- ✓ Three-dimensional stretch substrate
- ✓ Superior strength
- ✓ Weight bearing in approximately 20 minutes
- ✓ Excellent end lay down
- ✓ Reduces tucks or folds
- ✓ Superior lamination
- ✓ Easy open pouch
- ✓ 5-7 minute set time
- ✓ Available in ten popular colors

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)



# CAST REMOVAL SYSTEMS

## American Orthopedic™

Power & performance at the right price



A/O Cast Saw  
#0295-200

Specially designed for rugged use; perfect tool with power and endurance. Balance of weight and power helps reduce user fatigue. Internally lubricated for longer service. Switch is designed and placed so that it will not break if dropped. Hex drive provides 6 blade rotations extending blade life by 50% compared to similar competitive models. Length 12.5" Weight 3 lbs 13oz.

*The A/O Cast Saw is compatible with your choice of vacuum*



A/O Mobile Vacuum  
#0295-400

AMERICAN ORTHOPAEDIC™ Vacuum's quiet operation eases patient fears through reduction of noise. However, it is powerful enough to pick up even fine particles. M-PACT vacuums are the only vacuums that incorporate a unique filter cartridge thus creating a safer environment. Once the dust enters the filter it is trapped. Other vacuums use paper bags that are messy to clean. The remote on/off switch gives the ability to use the vacuum alone or with the saw.



A/O Portable Vacuum  
#0295-300

AMERICAN ORTHOPAEDIC™ Portable Vacuum's unique filtration system collects cast dust in a disposable filter cartridge- replacement is quick, easy and clean with no more messy bags. The convenient storage within the self-contained unit has space for cast spreader, scissors, replacement blades and the filter storage caps. The lightweight construction provides easy mobility. Furthermore, the portable vacuum is versatile as it attaches easily to most cast cutters.

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)



# CAST REMOVAL SYSTEMS

## Turbo Care

### The Gold Standard in Cast Removal Systems



TurboCare Cast Saw w/Power supply  
#TCCSS1

TURBOCARE® Cast Saw. Designed to reduce noise to a minimum results in less emotional trauma for patients and less hearing trauma for the clinician. The innovative motor plus saw design delivers high power in a smaller, lighter package. In addition, the unique design enhances the saw's balance and ergonomics for more user comfort. Plus automated multi-speed controls with "turbo-mode" for increased power.




TurboCare Vacuum  
#TCCVS1

TURBOCARE® VACUUM. The vacuum's quiet operation eases patient fears through reduction of noise. However, it is powerful enough to pick up even fine particles. M-PACT vacuums are the only vacuums that incorporate a unique filter cartridge thus creating a safer environment. Once the dust enters the filter it is trapped. Other vacuums use paper bags that are messy to clean. The remote on/off switch gives the ability to use the vacuum alone or with the saw.

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# CAST WEDGES



CastWedges aid in realignment and accommodate swelling. When used after an initial injury or surgery, the CastWedge permits uni- or bi-valving the cast, preventing extreme discomfort due to swelling. The CastWedge will maintain a precise opening in the cast. For realignment, CastWedge allows for exact adjustments to casted fractures using the opening wedge cast adjusting technique. The opening wedge cast/fracture technique avoids the inconvenience and cost of replacing the original cast. CastWedge is lightweight, inexpensive and can be used with plaster or fiberglass casts. It is made of high-impact, non-porous, radiolucent material that will not interfere with X-rays.



The CastWedge starter set consists of 100 wedges, 20 of each size (10mm, 15mm, 20mm, 25mm and 30mm). Replacement wedges are available in packages of 20 wedges.



Order #990-00

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# CERVICAL

## Aspen Collars

Aspen Cervical Collars are designed to restrict cervical spine flexion, extension and rotation to promote patient recovery. Two-piece polyethylene design enhances patient comfort and is easily adjusted with Velcro™ tabs. Anterior and posterior flex tabs on the front and back panels provide improved rigidity and durability. The large tracheal aperture provides functional access for the airway management and other essential procedures. Removable, foam laminated pads wick away moisture to ensure proper patient hygiene and comfort while preventing skin irritation. X-Ray, CT and MRI lucent. Innovative height adjustment technology on VISTA and VISTA TX offers six (6) collars in one. Collar Set includes extra set of replacement pads.



**Vista® Collar**



**Vista® TX Collar**



**Aspen® Collar**

*Order and sizing information on back*

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# CERVICAL



**Vista® Collar**

Part Number:	Description	Circ.	Size
79-83380	Vista Collar set	13"-19"	Adult
79-83370	Vista Collar	13"-19"	Adult
79-83371	Vista Collar repl pads		Adult
79-83377	Vista Collar back		Large



**Vista® TX Collar**

Part Number:	Description	Circ.	Size
79-83360	Vista TX Collar Set	13"-19"	Adult
79-83390	Vista TX Collar	13"-19"	Adult
79-83371	Vista Collar repl pads		Adult
79-83377	Vista Collar back		Large



**Aspen® Collar**

Part Number:	Description	Height/Age	Neck Circ.
79-83363	Adult Short-Collar set	2.25"	13"-22"
79-83365	Adult Regular- Collar set	3.0"	13"-22"
79-83367	Adult Tall- Collar set	3.75"	13"-22"
79-83368	Adult XTall- Collar set	4.5"	13"-22"
79-83283	Adult Short	2.25"	13"-22"
79-83285	Adult Regular	3"	13"-22"
79-83287	Adult Tall	3.75"	13"-22"
79-83288	Adult XTall	4.5"	13"-22"
79-83272	Infant Short	1-18 months	8"-12"
79-83273	Infant Regular	9-24 months	11"-13"
79-83275	Child Tot	1-3 years	11.5"-14"
79-83277	Child Short	2-5 years	12"-14.5"
79-83278	Child Regular	3-6 years	13"-16"
<b>Accessories</b>			
79-83295	Replacement Pads	Child	
79-83325	Replacement Pads	Adult	
79-83383	Back	S	
79-83387	Back	L	

**Southeast Medical Products, Inc**  
**5524 Commerce Dr**  
**Orlando FL 32839**

**407.851.5848 (ph) • 407.851.9012 (fax)**  
**800.330.0890 (ph) • 800.851.9012 (fax)**

# CERVICAL

## Miami J® Collar

Patented and proven, the Miami J® meets the rehabilitation needs of adults. Based on Clinical Studies and the input of leading medical institutions, this ground-breaking design provides the best immobilization to optimize patient rehabilitation, protects skin integrity during extended wear, and offers comfort features that enhance compliance and help speed recovery. All of our collars are latex-free and compatible with all imaging modalities.



**Miami Jr.®**

Pediatric c-spine care is distinct from adult care. The reason is basic — children are built differently. The Miami Jr.® from Jerome is the breakthrough collar designed exclusively for children 12 and under. The Miami Jr. line spans the growth cycle, is sized according to the Broselow-Luten color-coding system and meets the developing anatomical requirements unique to children.

When you snap on the Miami JTO Thoracic Extension, you make a great collar even better. Designed in conjunction with leading spine surgeons, orthotists and nurses, the JTO is used exclusively with all Miami J collars when even greater stability at the cervical junctions is indicated.



**Miami JTO®  
Thoracic Extension**

**Southeast Medical Products, Inc  
5524 Commerce Dr  
Orlando FL 32839**

**407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)**

# MIAMI COLLARS

## PRODUCT INFORMATION

Effective March 1, 2005



Size	Miami J <sup>®</sup> Collar Sets (collar + replacement pads)	Miami J <sup>®</sup> Collar	Miami J <sup>®</sup> Replacement Pads
<b>Specialty Sizes:</b>			
Super Short**	MJR-200S	MJ-200S	MJP-100
XS	MJR-250	MJ-250	MJP-250
Stout	MJR-200L	MJ-200L	MJP-200L
<b>Standard Sizes:</b>			
Small	MJR-300	MJ-300	MJP-100
Medium	MJR-400	MJ-400	MJP-100
Large	MJR-500	MJ-500	MJP-100

Size	Miami J <sup>®</sup> Front Panels	Miami J <sup>®</sup> Back Panels
<b>Specialty Sizes:</b>		
Super Short**	MJ-200S FRT	MJ-400 BK
XS	MJ-250 FRT	MJ-250 BK
Stout	MJ-200L FRT	MJ-200L BK
<b>Standard Sizes:</b>		
Small	MJ-300 FRT	MJ-300 BK
Medium	MJ-400 FRT	MJ-400 BK
Large	MJ-500 FRT	MJ-500 BK

Miami J<sup>®</sup> Front and Back Panels may be interchanged to accommodate patient anatomy, e.g., a Small Front may be used with a Large Back.

\*\*Super Short collars are packaged with Medium Back Panels



Size	Miami Jr. <sup>®</sup> Collar Sets (collar + replacement pads)	Miami Jr. <sup>®</sup> Collar	Miami Jr. <sup>®</sup> Replacement Pads
0 – 6 months*	MJR-P0	MJ-P0	MJP-P0
6 months – 2 years	MJR-P1	MJ-P1	MJP-P1
2 – 6 years	MJR-P2	MJ-P2	MJP-P2
6 – 12 years	MJR-P3	MJ-P3	MJP-P3

\* P0 collars are packaged with P1 Back Panels



Miami Occian Collar Back	ACB-100
--------------------------	---------

Southeast Medical Products, Inc  
 407.851.5848/800.330.0890 (ph)  
 407.851.9012/800.851.9012 (fax)  
 orders@southeastmedical.com

# CERVICAL

## Philadelphia Tracheotomy® Collar

Anterior and posterior rigid plastic reinforcement supports limit movement. Total cervical arch support maintains neutral alignment of the cervical spine.

- Water-resistant for bathing and aqua-therapy
  - Easily customized for individual requirements
  - Latex-free, non-toxic and hypoallergenic
- Plastazote® foam reduces skin irritation



	Neck Circ ↓	1.25"	1.75"	2.25"	3.25"	4.25"	5.25"
<b>Infant</b>	6-8"	PHP-TI					
<b>Pedi</b>	8-11"		PHP-TP				
<b>Xsmall</b>	Up to 10"			PHP-T2XS	PHP-T3XS	PHP-T4XS	PHP-T5XS
<b>Small</b>	10-13"			PHP-T2S	PHP-T3S	PHP-T4S	PHP-T5S
<b>Medium</b>	13-16"			PHP-T2M	PHP-T3M	PHP-T4M	PHP-T5M
<b>Large</b>	16-19"			PHP-T2L	PHP-T3L	PHP-T4L	PHP-T5L
<b>Xlarge</b>	19" +			PHP-T2XL	PHP-T3XL	PHP-T4XL	PHP-T5XL

Southeast Medical Products, Inc  
 5524 Commerce Dr  
 Orlando FL 32839  
 407.851.5848 (ph) • 407.851.9012 (fax)  
 800.330.0890 (ph) • 800.851.9012 (fax)



# C E R V I C A L

## Pneu-trac® Cervical Traction Collar

A comfortable, lightweight and easy-to-use alternative to traditional over-door traction, the Pneu-trac® provides cervical stretch through an inflatable collar. Clients can vary air pressure appropriate to the amount of cervical distraction required by using the hand-held inflation bulb with an air release valve. Made of semi-rigid polyethylene. Comes in three sizes.



- Pneumatic distraction control
- Patient adjustable
- Lightweight, semi-rigid polyethylene construction
- Two piece construction allows for varying height adjustments
- Velcro closure

Re-oder #	Size	Neck circ.
A26700	Small	12¼ - 14¼"
A26701	Medium	14½ - 17"
A26702	Large	17 - 20"

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# D M E - CRUTCHES/CANE/WALKERS



Lightweight aluminum crutches feature a double extruded center tube that provides additional strength to weight-bearing area. Push-pin adjustment makes it easy to adjust crutch height in 1" increments. Non-skid, jumbo size vinyl contoured tips provide excellent traction.

Item #	Size	Patient height
10400	Adult	5'2" – 5'10"
10401	Youth	4'6" – 5'2"
10402	Tall	5'10" – 6'6"
10416	Pediatric	4' – 4'6"



Adjustable aluminum cane features ergonomically designed handle with soft foam grip for comfort and security. Height adjusts from 30" to 39". Manufactured with sturdy, 1" diameter anodized extruded aluminum tubing.

Item #	Color	Adjustable Height
10307	Bronze	30" – 39"

### Quad Cane

Item #	Color	Base
10300	Chrome	Large
10301	Chrome	Small



Aluminum walker features single-button release to allow quick, one-hand folding. Cross brace is located near hand grips to allow for full stride. Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight. Includes rear glide caps, (item # 10107) and glide covers (item #10107C) allowing use on all surfaces.

Item #	Size	Adjustable Height
10202	Adult	32" – 39"
10203	Youth	25" – 32"

5" wheels, item#10109



**Southeast Medical Products, Inc**  
**5524 Commerce Dr**  
**Orlando FL 32839**

**407.851.5848 (ph) • 407.851.9012 (fax)**  
**800.330.0890 (ph) • 800.851.9012 (fax)**

# DME - FOREARM CRUTCHES



- Vinyl coated, ergonomically contoured arm cuffs molded for comfort and stability.
- Leg and forearm sections adjust independently for optimal sizing.
- Vinyl hand grips are comfortable and durable.
- Extra large tips for added stability.
- Ortho K handle is hard plastic.
- Limited Lifetime Warranty.

Item #	Description	Patient height	Handle height	Weight limit
10403	Adult	5' 0" – 6' 2"	28" – 37"	300 lbs
10404	Youth	4' 2" – 5' 2"	21.5" – 30.5"	300 lbs
10405	Tall	5' 10" – 6' 6"	33" – 42"	300 lbs
10407B	Child, Blue	3' 2" – 4' 6"	14" – 21"	250 lbs
10407R	Child, Red	3' 2" – 4' 6"	14" – 21"	250 lbs

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# DORSAL CARPAL TUNNEL

The RCAI Dorsal Carpal Tunnel Splint is made of a lightweight, durable and heat moldable thermoplastic for patient customization. The narrow palmar roll supports the hand while fingers remain free to grasp.



The Dorsal Carpal Tunnel Splint limits wrist flexion and extension during repetitive hand motion. It can also be used as a wrist extension splint to support a weak or flaccid wrist.

The breathable foam liner wicks moisture away from the skin.

<b>Size</b>	<b>Measure width of palm</b>	<b>Model # Right</b>	<b>Model # Left</b>
Ex Small	2" - 2 1/2" (5.1 - 6.4 cm)	44DCT-ES-R	44DCT-ES-L
Small	2 1/2" - 3" (6.4 - 7.6 cm)	44DCT-S-R	44DCT-S-L
Medium	3" - 3 1/2" (7.6 - 8.9 cm)	44DCT-M-R	44DCT-M-L
Large	3 1/2" - 4" (8.9 - 10.2 cm)	44DCT-L-R	44DCT-L-L
X-Large	4" - 4 1/2" (10.2 - 11.4 cm)	44DCT-XL-R	44DCT-XL-L

Southeast Medical Products, Inc  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# EDEMA GLOVES

Soft, comfortable compression gloves for edema patients.  
Available in 3/4" or full finger length. Color: Champagne.



Part	Description		
	Length	Size	Hand
902XSL	7"	X-Small	Left
902XSR	7"	X-Small	Right
902SL	8"	Small	Left
902SR	8"	Small	Right
902ML	9"	Medium	Left
902MR	9"	Medium	Right
902LL	10"	Large	Left
902LR	10"	Large	Right

Part	Description		
	Length	Size	Hand
903XSL	7"	X-Small	Left
903XSR	7"	X-Small	Right
903SL	8"	Small	Left
903SR	8"	Small	Right
903ML	9"	Medium	Left
903MR	9"	Medium	Right
903LL	10"	Large	Left
903LR	10"	Large	Right



Part	Description		
	Length	Size	Hand
904SL	8"	Small	Left
904SR	8"	Small	Right
904ML	9"	Medium	Left
904MR	9"	Medium	Right
904LL	10"	Large	Left
904LR	10"	Large	Right

Specify RIGHT or LEFT. Measure around the hand, palm flat, around the knuckles excluding the thumb. Order larger size if between sizes.

**Southeast Medical Products, Inc**  
**5524 Commerce Dr**  
**Orlando FL 32839**  
**407.851.5848 (ph) • 407.851.9012 (fax)**  
**800.330.0890 (ph) • 800.851.9012 (fax)**

# ELECTROTHERAPY

## Intelect® Portable Electrotherapy

*For symptomatic relief and management of chronic pain and acute post-surgical or post-trauma pain.*



- Dual, independent channels
- Adjustable parameters
- 5 stimulation modes (Burst, Normal, MRW, SD, Bi-pulse)
- Patient compliance timer
- Digital LCD display

**#77712 Intelect® TENS Digital**



- Dual, independent channels
- Adjustable parameters
- 3 stimulation modes (Burst, Normal, Modulation)
- Flip cover to protect top controls

**#77600 Intelect® TENS Standard**



- Dual, independent channels
- Adjustable parameters
- 3 stimulation modes (Burst, Normal, Modulation)

**#77710 Intelect® TENS Basic**

**DURA-STICK® II**  
SELF-ADHESIVE ELECTRODES



*Our most popular and best selling electrode, DURA-STICK® II electrodes provide comfortable and effective stimulation through the unique silver printed conductor.*

Flexible, yet sturdy construction offers increased durability and ease of removal. Soft and conformable lead wire connectors offer superior grip to the lead wire pin.

*All devices come with lead wires, 9-volt battery, DURA-STICK® II Self-Adhesive Electrodes and a convenient carrying case.*

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# FINGER SPLINTS



STAX SPLINTS are designed to support distal joint of finger in extension while permitting unrestricted movement of proximal interphalangeal joint. Available in starter kit and individual sizes (sold each).



Digital splints in 11 pre-cut sizes meet most splinting needs. Available in starter kit and individual sizes (sold by the dozen)



Ceramic finish malleable aluminum strips with radius ends that can be cut to length needed. Available in blue, soft, open cell foam and white, firm, non-absorbent, closed cell foam. Sold by the dozen.



SPECIALTY KIT. An assortment of 45 popular padded specialty splints in a plastic organizer.



COT Aluminum splint folds evenly over the end of finger to provide stabilization and protection.



PROTECTOR Aluminum splint provides end protection to injured fingers.



SPOON Aluminum splint curves upward at tip to protect end of finger.



FROG Aluminum splint folds to hold finger in proper alignment.



BASEBALL Aluminum splint folds over middle of finger and easily adjusts and molds to finger.

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)



# FINGER SPLINTS



## 3pp™ ThumSling™

Gently hugs the thumb and applies light compression to reduce grinding and pain at the CMC (basal) joint of the thumb.



### Features

- Wrap around sling repositions the CMC joint at the base of the thumb so you can perform daily activities with less pain or fear of dropping things.
- Easy to apply, this soft, flexible splint allows you to grip a golf club, hold a rake or grip a pen securely and comfortably.
- Lightweight and lined in soft, breathable foam, the ThumSling is not hot or bulky. Machine washable.
- Light Control

### Sizing Information

Instructions: Lay hand flat and measure around the MP joints (large knuckles).

Item Number	Size	Fits
P2006	Small/Medium	6" to 7 1/4"
P2006	Medium/Large	7 1/4" to 9 1/4"

Choose right or left hand.

## 3pp™ Buddy Loop™

Eliminates messy, sticky tape and prickly straps! Buddy Loops are the number one choice of physicians and athletes to protect finger injuries at rest and during play.



### Features

- Fast "loop and wrap" design makes it easy to hold fingers together securely.
- Soft, foam-lined material provides a light cushion between your fingers for comfort and control.
- Use one or two Buddy Loops as needed for control. For larger hands or longer fingers, we recommend a wider one-inch strap at the base of the fingers used in combination with a half-inch strap near the fingertips.
- Machine washable.
- Light Control

### Sizing Information


Instructions: 5" long loops fit most fingers. Available in three widths - 1/2" wide, 3/4" wide and 1" wide. Sold in packages of 5, 25, 50 or 100.

Item Number	Size
P1005	1/2" Buddy Loop
P1007	1" Buddy Loop
P1003	3/4" Buddy Loop

Southeast Medical Products, Inc  
5524 Commerce Dr  
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# FOOT PRODUCTS



## Rolyan® External Fixator Foot Support

A dynamic dorsiflexion assist for tibial fractures. Helps prevent foot drop and protects heel. Compatible with external fixation systems: including Hoffman, Kroner and Fisher. Elastic straps give dorsiflexion support and plantar flexion exercise. Can be positioned for inversion or eversion. Base washable.

### External Fixator Replacement

Includes replacement foam, two elastic straps with snap catch, and one foam strap pad.

Item#	Description
A4981	Rolyan® External Fixator Foot Support
A4983	Rolyan® External Fixator Foot Support External Fixator Replacement

Southeast Medical Products, Inc  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# FOOT PRODUCTS / SPENCO



## PolySorb Cross Trainer Insoles #38-034-xx

*For all sports or athletic activity*

For aerobics, court sports, running, weight training, racquet sports and walking. Special cushioning in both the heel strike zone and the forefoot area for outstanding performance. Helps turn any athletic shoe into a cross training shoe.



## Slip-In Insoles #40-212-xx

*Long lasting heel-to-toe comfort*

Years of proven use and high quality, synonymous with the name SPENCO®, our slip-in insoles provide long lasting heel-to-toe comfort. Washable, trim to fit with scissors.



## 3/4 Length Orthotic Arch Supports #43-158-xx

Spenco® Orthotic Arch Supports are specifically designed to help control the position of the feet. They are shaped to help stabilize and balance the foot and to provide the extra comfort and support needed for optimum performance. This unique cushioning system helps absorb shock, reduce friction and improve overall foot comfort. Also available in full length.



## Size Chart

xx	00	01	02	03	04	05	06
Women's	3/4	5/6	7/8	9/10	11/12		
Mens's			6/7	8/9	10/11	12/13	14/15

*All Spenco products are available from Southeast Medical Products.*

**Southeast Medical Products, Inc**  
**5524 Commerce Dr**  
**Orlando FL 32839**  
**407.851.5848 (ph) • 407.851.9012 (fax)**  
**800.330.0890 (ph) • 800.851.9012 (fax)**

# HAND SPLINTS

## Galveston Metacarpal Brace



- **ADVANTAGES:** Utilizes the concept of three point fixation, a better degree of correction than with the ulnar gutter splint, fully adjustable pads to customize the splint to various hand sizes and fracture types, and an adjustable strap that allows the maintenance of swelling reduction. It is radiotransparent, and allows unrestricted motion of IP, MP, and wrist joints.
- **INDICATIONS:** Can be applied to fractures of the 2nd through the 5th metacarpal on either hand by simply adjusting the crossbar and pads.
- **PRECAUTIONS:** The two-toned dorsal pad has a soft layer in contact with the skin, and a stiffer backing pad. As long as the soft pad is not fully compressed, pressure under the pad is below that shown to cause necrosis (the soft pad is only meant to be an indicator). As with other non-invasive methods of treatment, this splint does not prevent shortening, The Galveston Metacarpal Brace, as with any other functional bracing system, requires patient compliance for optimal results.
- **APPLICATION CONSIDERATIONS:** The brace should be applied following fracture reduction to maintain alignment. The patient should be allowed to adjust the brace to be snug, but not uncomfortable.

### Sizing

Measure the Diameter of the Palm.



Order #	Description	Inches	Centimeters
GMB101S	Small	2.25"-3"	6-8cm
GMB101M	Medium	3"-3.5"	8-9cm
GMB101L	Large	3.5"-4"	9-10cm
GMB101XL	Extra Large	Over 4"	Over 10cm

A spare set of pads is included with each brace.  
Extra pads & bars are available at an additional cost.

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# HEEL CUPS & PADS



## Silicone Heel Cups

Molded medical-grade soft silicone to help reduce shock and vibration loads to joints, to address plantar fasciitis, heel spurs and achilles tendonitis. Softer center in heel region to help reduce strain to area sensitive to pain. Sold in pairs.

Number:	Male	Female	Size
79-81102	Children/youth	Children/youth	XS
79-81103	5-1/2 - 9-1/2	5 - 9	S/M
79-81107	10+	9-1/2+	L/XL



**Aetrex Heel Spur Cradles** were doctor designed to comfort heel pain. The cradle absorbs shock and cushions painful heels for immediate relief. Anti-Shox Spur Heel Cradle provides relief from heel spurs and plantar fasciitis.

Number:	Size
5075S	SML
5075M	MED
5075L	LRG



## Bone Spur Pad

To reduce pain in the heel. Outer ring absorbs shock while inner circle provides an area to gently accommodate injury by containing soft tissue around heel. Relieves pressure on bone spur or in sensitive area of heel. Can help relieve lower back fatigue, heel spurs, plantar fasciitis and general heel discomfort.

Number:	Male	Female	Size
71301	4 - 6	5 - 7	SML
71302	6 1/2 - 10 1/2	7 1/2 - 11 1/2	MED
71303	11 +	12 +	LRG



## Tuli's® Classic Heel Cups

Protect your heel in biomechanically-designed cradles that provide relief from heel pain.

Number:	Size
T10010	REG
T10020	LRG



**Tuli's® Heavy Duty Heel Cups** Equip your heel with the ultimate shock absorption required for high impact work or play.

Number:	Size
T10211	SML
T10212	MED
T10214	LRG



## TuliGEL™ Heel Cups

Tuli's classic support in a lighter, softer but resilient gel cup.

Number:	Size
T10215	REG
T10216	LRG

**Southeast Medical Products, Inc**  
**5524 Commerce Dr**  
**Orlando FL 32839**

**407.851.5848 (ph) • 407.851.9012 (fax)**  
**800.330.0890 (ph) • 800.851.9012 (fax)**

# L / S SUPPORTS - Warm/Form style



## **79-8250x**

Low profile 9" elastic support with double side pull compression straps. Neoprene posterior panel accommodates ComfortFORM™ Moldable Insert #79-89120. Foam lumbar pad is included with 79-8250x style. No stays. Ideal for providing support and compression for strains & sprains of the low back and abdominal area.



## **79-8925x**

9" neoprene/elastic construction. Low-profile design. Contact closure. Accommodates COMFORT-FORM Moldable Insert #79-89120.



## **79-8924x**

Foam/pile, with 9" panel elastic construction and double side-pull elastic compression straps. Two sizes. Accommodates ComfortFORM™ moldable insert #79-89120. Ideal for providing support and compression for strains, sprains of the low back and abdominal area.



## **79-8249x**

9" all-neoprene construction with double side-pull compression straps for additional support. 2-ply posterior pocket accommodates ComfortFORM™ Moldable Insert #79-89120. Ideal for providing support and compression for strains & sprains of the low back and abdominal area.



The ComfortFORM™ Moldable Insert can be heated with a heat gun or at 160° in a hot water pan, or hydrocollator and then formed against patient's back for exact contour and improved support. To be used with a compatible lumbosacral support to provide rigid support to the lower back.

**Southeast Medical Products, Inc**  
**5524 Commerce Dr**  
**Orlando FL 32839**

**407.851.5848 (ph) • 407.851.9012 (fax)**  
**800.330.0890 (ph) • 800.851.9012 (fax)**




# L / S S U P P O R T S

## Saunders Posture S'port®



**The Posture S'port** is unique because it anchors at the waist and not only discourages stooped shoulders, but assists in proper spinal alignment. The Posture S'port can be worn directly on the skin or over a t-shirt.



### Features:

- Criss-crossed suspender like straps gently pull the shoulders back into the neutral position to decrease back and shoulder pain
- Soft cotton under arm pads provide a comfortable fit that won't pinch or bind
- Can be worn directly on skin or over clothing
- Adjustable elastic waistband that anchors the entire support to the waist

Small	22" - 30"
Medium	30" - 42"
Large	42" - 54"

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)



# MASSAGE THERAPY

## Adapta® MC-100 Portable Massage Chair



Portability is a prerequisite for any successful massage therapist.  
Available in four colors: Royal Blue (#92019), Wedgewood (#92022),  
Maroon (#92020) and Mushroom (#92021).

Lightweight - 16 lbs. (7 kg)

UL listed at 300 lbs. (135 kg) working weight capacity

Durable and affordable

**Shown with OPTIONAL Deluxe Carry Case, #93001.**

Southeast Medical Products, Inc  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# PLANTAR FASCIITIS SPLINTS



**procare™**

**Dorsiwedge Night Splint** - Soft padded, contoured posterior shell provides passive dorsi-flexion to address plantar fasciitis and Achilles tendonitis. Fits right or left foot. Shipped with toe wedge for additional stretch on plantar fascia.



**DARCO.**



The **Darco Night Splint** is designed to alleviate heel pain associated with plantar fasciitis and achilles tendonitis.

Special Features:

- Terry cloth lining is all natural, comfortable and moisture absorbent
- Non-slip tread on the sole
- Pneumatic dorsi-flexion air bladder provides infinite angle adjustment
- Maintains proper flexation when worn properly

**BIRD & CRONIN INC.**



## **Baker™ Plantar Fasciitis Night Splint**

- *Designed to help minimize the pain and discomfort of plantar fasciitis*
- *Allows the patient to apply a consistent 90° angle stretch to the plantar fascia by holding the foot in gentle dorsiflexion*
- *Soft padded cover provides greater patient comfort and protection against irritation to the skin and breakdown*
- *Removable foam wedge insert is included for additional 5" stretch on the plantar fascia*
- *Fits either left or right foot*
- *Latex Free*

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# PLANTAR FASCIITIS SPLINTS



## P.F.S.™ Plantar Fasciitis Stretch Splint

- *Designed to provide relief from pain and discomfort of plantar fasciitis and Achilles tendonitis*
- *Allows the patient to apply consistent stretch to the plantar fascia by holding the foot in gentle dorsiflexion*
- *Lightweight design for patient comfort while sleeping*
- *Slip resistant sole for greater patient safety*
- *Fits either left or right foot*
- *Latex Free*



## Formfit Night Splint

- Adjustable flexion straps provide a gentle stretching of the plantar fascia and achilles tendon. Adjusts from 90 to 10 dorsi-flexion
- Patented 3-dimensional molding creates superior comfort and fit with padding where you need it
- Soft, lycra lined foam interior is plush, washable and breathable to eliminated heat build-up
- Lycra lined foam padded straps maintain comfort across the flexor tendons
- User-friendly center-release buckles simplify application
- Lightweight, low profile shell is sturdy and breathable



**Southeast Medical Products, Inc**  
**5524 Commerce Dr**  
**Orlando FL 32839**

**407.851.5848 (ph) • 407.851.9012 (fax)**  
**800.330.0890 (ph) • 800.851.9012 (fax)**

# PLANTAR FASCIITIS SPLINTS



## AirForm Dorsal Night Splint

### Performance Benefits:

- The Airform Night Splint provides an effective treatment for plantar fasciitis, achilles tendonitis, drop foot and post-static pain
- Its lightweight, low profile design makes it more comfortable. Better comfort means better compliance
- Slips on and adjusts easily
- Fits left or right foot
- Comes in 2 sizes

### Indications:

- Plantar Fasciitis
- Achilles Tendonitis
- Drop Foot
- Post-Static Pain



## The Strassburg Sock

For treatment of Plantar Fasciitis, heel spurs, Shin Splints and Achillies tendonitis.

### Features:



- It is very easy to use.
- It does not lock the ankle joint into any one position thus allowing some movement between the various joints preventing possible joint compression pain.
- It can be safely worn for several hours without adverse effect when applied properly.
- It can be worn to bed, provides no significant restrictions and will not interrupt normal sleep if worn as directed.
- It is an effective device providing low cost relief of heel spur symptoms.
- It can reduce/eliminate the pain experienced with the first morning steps.
- It can, over a period of time, strengthen the arch of the foot.
- It can supplement and not interfere with the various traditional doctor's treatment.

**Southeast Medical Products, Inc**  
**5524 Commerce Dr**  
**Orlando FL 32839**

**407.851.5848 (ph) • 407.851.9012 (fax)**  
**800.330.0890 (ph) • 800.851.9012 (fax)**

# REHABILITATION PRODUCTS

## 3x10=0 Full Extension Knee Device



The 3x10=0 Full Extension Knee Device is for the patient who lacks passive extension of the knee following ACL reconstruction, total knee replacement, knee surgery, and/or trauma. Upon first observation of knee extension lag, this prophylactic appliance allows patients to regain full extension through a safe and effective rehabilitation routine.

With the thin end under the thigh, the ankle is supported by the riser at the distal end of the foam cushion. The straps are drawn around the leg and the foam cushion. The straps are then tightened gradually to move the knee toward full extension.

The 3x10=0 Full Extension Knee Device is easier to use than the alternative day and night devices, therefore patient compliance is high. It can also be advantageous in overcoming long standing extension problems by simply lengthening the application time.

**Southeast Medical Products, Inc**  
**5524 Commerce Dr**  
**Orlando FL 32839**  
**407.851.5848 (ph) • 407.851.9012 (fax)**  
**800.330.0890 (ph) • 800.851.9012 (fax)**

# RESOLVE® HALO

Ossur's ReSolve Halo System is the industry's only MR-Safe halo system. ReSolve combines a glass composite ring with Ossur's proprietary ceramic-tipped skull pins.

It is the only halo clinically proven to eliminate the risk of dangerous scalp burns resulting from halo pin heating phenomena.

ReSolve isolates the patient from potential current induction in the [MRI](#). Not only does the Resolve system ensure MR-safety but clearer and more high-resolution brain images are now possible.



## Features and Benefits

- The industry's only MR-Safe Halo System
- The only halo clinically proven to eliminate the risk of dangerous scalp burns resulting from halo pin heating phenomena
- Glass composite ring and ceramic-tipped pins are non-conductive and induce no dangerous current flow
- MR images are clear and artifact-free

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# RESOLVE® HALO

Part Number	Description
Resolve Nonconductive Halo Rings	
505300C	Resolve Closed Back Halo • 17-20" (43-51cm), Small
505400C	Resolve Closed Back Halo • 20-23" (51-59cm), Standard
505500C	Resolve Closed Back Halo • 23-26" (59-66cm), Large
510300C	Resolve Closed Back Halo, Ring Only • 17-20" (43-51cm), Small
510400C	Resolve Closed Back Halo, Ring Only • 20-23" (51-59cm), Standard
510500C	Resolve Closed Back Halo, Ring Only • 23-26" (59-66cm), Large
505300D	Resolve Open Back Halo • 17-20" (43-51cm), Small
505400D	Resolve Open Back Halo • 20-23" (51-59cm), Standard
505500D	Resolve Open Back Halo • 23-26" (59-66cm), Large
510300D	Resolve Open Back Halo, Ring Only • 17-20" (43-51cm), Small
510400D	Resolve Open Back Halo, Ring Only • 20-23" (51-59cm), Standard
510500D	Resolve Open Back Halo, Ring Only • 23-26" (59-66cm), Large
Halo Vests	
545200C	Halo Vest Assembly • 20-26" (50-66cm), X-Small
545300C	Halo Vest Assembly • 26-33" (66-84cm), Small
545350C	Halo Vest Assembly • 32-36" (81-92cm), Small/Medium
545400C	Halo Vest Assembly • 33-40" (84-102cm), Medium
545450C	Halo Vest Assembly • 38-42" (96-107cm), Medium/Large
545500C	Halo Vest Assembly • 40-46" (102-117cm), Large
545600C	Halo Vest Assembly • 46-52" (117-132cm), X-Large
Lil' Angel Pediatric Vests	
545100K	Lil' Angel Vest, P2 w/ Kodel • P2 (2-6 Years)
545100L	Lil' Angel Vest, P2 w/ Lambswool • P2 (2-6 Years)
545100S	Lil' Angel Vest, P2 w/ Sorbatex • P2 (2-6 Years)
545150K	Lil' Angel Vest, P3 w/ Kodel • P3 (6-12 Years)
545150L	Lil' Angel Vest, P3 w/ Lambswool • P3 (6-12 Years)
545150S	Lil' Angel Vest, P3 w/ Sorbatex • P3 (6-12 Years)

**Southeast Medical Products, Inc**  
**5524 Commerce Dr**  
**Orlando FL 32839**  
**407.851.5848 (ph) • 407.851.9012 (fax)**  
**800.330.0890 (ph) • 800.851.9012 (fax)**



# SHOES, CAST SHOES



**Available in 11 sizes!**



**#882 Cast Shoe** - Our most popular style upper, made of strong canvas duck. Straps of canvas duck have hook and loop closures for a more secure fit. The rocker sole aids in the mobility of the patient and the open toe and heel allows for more versatility in sizing.

Order#	Size	sole length
882-98	3XS	4 7/8"
882-INF	2XS	5 1/2"
882-000	Todlr	6 3/8"
882-00	XSM	7 1/4"
882-005	XSM/SML	8"
882-01	SML	9 1/2"
882-015	SML/MED	10 1/2"
882-02	MED	11 1/4"
882-025	MED/LRG	12 1/4"
882-03	LRG	12 3/4"
882-04	XLG	13 1/2"

**The Darco SlimLine Cast Boot** has taken another step forward in providing cast protection in today's active world. Excellent for use with fiberglass casts or bulky compression bandages because the stretchable upper material conforms for a better fit.

Special Features:



- Square toe design acts as a protective bumper and gives the patient more room to accommodate larger casts and bandages.
- Rocker sole allows for smooth roll through the gait pattern.
- Removable EVA sole provides 10mm of padding to cushion the foot, adding to patient comfort.
- Strapless forefoot closure is easily adjusted to accommodate cast/bandage size and can be switched between left/right closure for ease of re-application.
- Elastic center strap secures neck of the ankle and seats the cast in the shoe.
- Higher ankle strap adds stability

Item No.	DESCRIPTION	SIZE	COLOR	WOMEN SIZE	MEN SIZE
SLQ0B	SlimLine Cast Boot	X-Small	Black	3 - 5	
SLQ1B	SlimLine Cast Boot	Small	Black	5.5 - 8	
SLQ2B	SlimLine Cast Boot	Medium	Black	8.5 - 10	6 - 8
SLQ3B	SlimLine Cast Boot	Large	Black	10.5 - 13	8.5 - 11
SLQ4B	SlimLine Cast Boot	X-Large	Black	13.5 +	11.5 -14

**Southeast Medical Products, Inc**  
**5524 Commerce Dr**  
**Orlando FL 32839**  
**407.851.5848 (ph) • 407.851.9012 (fax)**  
**800.330.0890 (ph) • 800.851.9012 (fax)**

# SHOES, POST-OP



**Darco's Med-Surg Shoe's** strapless forefoot closure eliminates buckle pressure on the forefoot and can expand to accommodate dressings up to 17" in circumference. Available accessories; removable peg insole and toe cover.



ITEM No. (WOMEN)	ITEM No. (MEN)	DESCRIPTION	SIZE	COLOR	WOMEN SIZE	MEN SIZE
MQPB		Med-Surg Shoe	Pediatric	Black		
MQW1B	MQM1B	Med-Surg Shoe	Small	Black	4 - 6	6 - 8
MQW2B	MQM2B	Med-Surg Shoe	Medium	Black	6.5 - 8	8.5 - 10
MQW3B	MQM3B	Med-Surg Shoe	Large	Black	8.5 - 10	10.5 - 12
	MQM4B	Med-Surg Shoe	X-Large	Black		12.5 - 14



## ProCare Post-Op Shoe

Constructed with E.V.A./rubber, no-flex rocker design sole and contoured heel. Low-profile foam/nylon mesh upper has contact closure. Protective tongue provides maximum comfort.



## ProCare Med/Surg Shoe

Designed for a variety of post-op and post-trauma applications. Foam mesh upper with semi-flexible skid resistant sole and contoured heel. Contact closure straps. Fits right or left foot.

Part Number:		Shoe Size	Size
79-90183	Male Loop/Lock Closure	7 - 9	S
79-90185	Male Loop/Lock Closure	9 - 11	M
79-90187	Male Loop/Lock Closure	11 - 13	L
79-90188	Male Loop/Lock Closure	13+	XL
79-90193	Female Loop/Lock Closure	4 - 6	S
	Female Loop/Lock		

Part Number:	Shoe Size	Size
79-81142	2 - 4 Female	XS
79-81143	4 - 6 Female	S
79-81145	6 - 8 Female	M
79-81147	8 - 10 Female	L
79-81132	5 - 7 Male	XS
79-81133	7 - 9 Male	S
79-81135	9 - 11 Male	M
79-81137	11 - 13 Male	L
79-81138	13+ Male	XL



**#442 Post-Op Shoe** – The Super Shoe II has the proper inner sole platform and cushioned rocker sole for natural gait. Durable upper with a soft white inner lining. Molded heel counter provides a smooth and stable contour around the heel. Unique removable tongue for optimum placement and comfort. Loop lock loop-thru closure for easy yet secure fit.



Order#	Size
399V	2XS
442-00	XSM
442-01	WOMENS SML
442-02	WOMENS MED
442-03	WOMENS LRG
442-11	MENS SML
442-12	MENS MED
442-13	MENS LRG
442-14	MENS XLG

**Southeast Medical Products, Inc**  
**5524 Commerce Dr**  
**Orlando FL 32839**

**407.851.5848 (ph) • 407.851.9012 (fax)**  
**800.330.0890 (ph) • 800.851.9012 (fax)**

# SPECIALTY SHOES , *Pressure Relief*



**The Darco Peg-Assist System** is built around a multi-purpose removable peg insole that effectively off-loads the plantar aspect of the foot.

Special Features:

- Pegs can be removed to off-load pressure from any area of the plantar surface of the foot.
- Constructed of an 18mm combination of Plastazote®, Poron®, and EVA to provide the perfect combination of shock absorption, off-loading and moldability.
- Peg Support Board adheres to black side of insole to keep pegs from collapsing and keep the off-loaded area intact.

**The OrthoWedge Shoe** has been clinically proven to effectively reduce weight bearing pressure from the forefoot to promote proper healing.

Special Features:

- Square toe design acts as a bumper to provide additional protection when K-wires are present and provides better left/right fit.
- Promotes healing by reducing weight from the forefoot and may be used for any condition from the metatarsal heads distally in which it is desirable to reduce body weight such as ulcers, infections, trauma and surgery.
- 15 degree dorsi-angle shifts weight back to the midfoot and heel.
- Ankle strap eliminates heel slippage and removable forefoot closure eliminates buckle pressure.
- Wide base for extra stability.
- 10mm removable insole fully customizable with over twice the padding.



Also available in HeelWedge version.



**The IPOS Forefoot Relief Orthosis** relieves the forefoot by avoiding any floor contact of the anterior aspect of the foot through a 10° dorsiflexed inclination of the heel and a shortened weight-bearing sole.

- 10° dorsiflexed inclination of the heel.
- Weight-bearing sole ends in the mid foot area.
- Beveled heel for better heel to toe transfer.
- Three closure straps for optimal fit.
- Inner surface completely padded and without seams.
- Non-slip sole for outside use.



Also available in heel relief version.

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# SHOULDER PRODUCTS


## Saunders Sully Shoulder Support

Traditional shoulder braces have forced you to sacrifice either stability or function. Now you can have both!


On the outside elastic straps attach with hook to the neoprene at any point, in any direction, and with any amount of force you choose. You can functionally stabilize, assist or restrict movement according to the specific needs of each athlete.

On the inside the perforated, breathable, rubber-like neoprene grips the skin so the support moves with the athlete. It's as if the stabilization straps are attached directly to the body.


The Sully comes complete with a highly adaptable, bifurcated neoprene strap, and two additional 4" neoprene straps. Worn by professional and amateur athletes in all sports! Measure biceps in a non-flexed position.



POSTERIOR INSTABILITIES



MUSCLE STRAINS (PECTORIAL ETC.)



ROTATOR CUFF DECELERATION



ANTERIOR INSTABILITIES



INFERIOR INSTABILITIES



MULTI-DIRECTIONAL INSTABILITIES

Southeast Medical Products, Inc  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# SPINAL HYPEREXTENSION

## ACE™ Brace

Three reasons the ACE Brace is simply the most comfortable, cosmetic hyperextension orthosis on the planet!

1. One size fits all backstrap
2. No sternal pressure - no throat impingement
3. Front closure - eases donning by arthritic hands



### ACE-ST

Hinged sternal assembly, hinged pubic assembly, universal backstrap and pad, front closure. Ideal for younger, stronger patient who can tolerate direct sternal pressure and is the most economical of our ACE Braces.



### ACE-PV

Pectoral pads on sponge rubber disc with velcro® fastener, hinged pubic assembly, universal backstrap and pad, front closure. Ideal for geriatric osteoporotic patient with prominent bony sternum, where orthotist desires to allow patient to shift pectoral pad position for comfort.



### ACE-300

A full featured universally sized hyperextension orthosis with one-size-fits-all back strap with innovative Y-tab and swivels on both sternal and pubic plates. Aluminum plates for rapid, precise fit.



### ACE-PS

Pectoral pads with adjustable silicone swivels, hinged pubic assembly, universal backstrap and pad, front closure. Ideal for geriatric osteoporotic patient with prominent bony sternum, where orthotist desires pectoral pads to remain exactly as fitted and is the most cosmetic of our ACE Braces.

Southeast Medical Products, Inc  
5524 Commerce Dr  
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)



# SPINAL HYPEREXTENSION

## C.A.S.H. Cruciform Anterior Spinal Hyperextension

- Pre-fabricated hyperextension orthosis
- Wide straps
- Fully adjustable height & width
- Available with horizontal or vertical posterior pad
- Cruciform does not make patient contact when worn



Sternal Pad Model

- #CHO100 Fixed sternal and pubic pads
- #CHO300 Hinged sternal and pubic pads
- #CHO400 Hinged sternal pad/fixed pubic pad



Pectoral Extension Model

- #CHO100P Fixed sternal and pubic pads
- #CHO300P Hinged sternal and pubic pads
- #CHO400P Hinged sternal pad/fixed pubic pad

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# SPLINT ROLL

## OCL® Splint Roll®-2 Plaster Splinting System



The OCL® Splint Roll®-2 Custom Splinting System is an improved plaster splinting product. Greater strength than other plaster splints. The DuoSet™ plaster gives faster initial set on outer layer. The foam and flannel covering is comfortable for the patient. The open cell foam prevents splint from migrating and improves air exchange for enhanced patient comfort. The flannel reduces moisture retention; ensures faster set times. Non-woven layer allows plaster crystals to lock in the dimpled areas and form additional reinforcement without added weight. Thermal barrier also doubles as a filtering system to prevent plaster from embedding in foam layer. Continuous roll configuration means less waste; cut what you need. Sizing is easy with different widths available. Surged sides provide added patient comfort. Innovative design facilitates custom splinting applications.

Order#	Size
SR215	2"x20'
SR315	3"x20'
SR415	4"x20'
SR515	5"x20'
SR615	6"x20'

### Accessories



Disposable buckets make clean-up easier. Economical. Perfect for immersing products; plaster splint roll, plaster bandage and splints, casting tape, synthetic splint roll.



Templates are an aid for preparing custom splints: simply select pattern, trace, prepare splint, cut and apply to patient.

Scissors are 9" razor-edge shears to allow for easy cutting of OCL® PolyLITE® Synthetic Splinting System.

OCL® Cutting Knife comes with 3 steel blades stored inside for cutting plaster splinting roll.



Plaster Splint Roll® Dispenser  
- 5 compartment



M-PACT® Elastic Bandage High quality elastic retains its shape. Double length bandages available in 4" and 6" widths.

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)



# WALKING BOOTS

**MaxTrax walkers** are designed to increase stability following trauma or post-operative procedures. Designed with a wider foot bed for comfort and stability, this low profile boot features a rocker bottom which helps promote natural gait and reduce plantar pressures. Cushioned inner and outer sole provides shock absorption to help improve patient comfort during ambulation. Malleable uprights breakdown if needed for custom fit.

## MaxTrax™ Walker



## MaxTrax™ Air Walker



MaxTrax Tall	MaxTrax Short	Size	Men	Women	MaxTrax Air Short	MaxTrax Air Tall
79-95322	79-95342	Xsmall	NA	Up to 4	79-95422	79-95412
79-95323	79-95343	Small	up to 5	4.5-6	79-95423	79-95413
79-95325	79-95345	Medium	5.5-10	6.5-11	79-95425	79-95415
79-95327	79-95347	Large	10.5-13.5	11.5-14.5	79-95427	79-95417
79-95328	79-95348	Xlarge	14-17	NA	79-95428	79-95418

Southeast Medical Products, Inc  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# WALKING BOOTS

## Wee Walker Boot

The **Bledsoe Wee Walker Pediatric Boot** is the first off-the-shelf low profile pediatric boot that offers superior comfort along with unsurpassed support. Only two sizes are offered thus reducing inventory cost. These sizes will fit left or right, ages 18 months to 5 years old.



### Features include:

- Promotes faster healing through strong support with maximum mobility.
- Designed for kids on the go. It's lighter and easier than conventional casts.
- Unique rocker design simulates a near normal gait.
- Keeps kids comfortable with breathable, washable foam wrap.
- Kid-tough performance from durable powder coated aluminum

Children's Shoe Size	Foot Length	Leg Length	Wee Walker Size	Part Number
up to 8	up to 6"	up to 9"	Small - Medium	AL185003
8-12.5	6"-7.5"	9"-12"	Large - XLarge	AL185007

### Indications

This device is indicated for post-injury or post-operative application to grade 1, 2, or 3 ankle sprains, and stable or internally fixed fractures of the foot or ankle. Also indicated for treatment of stable or internally fixed fractures or fractures with bridging callous of the middle to distal tibia or fibula.

### Contraindications

This device is contraindicated for unstable fractures and proximal tibial or fibular fractures.

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839  
407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)

# WRIST ORTHOSES (WH & WHT)

## Wrist Hand Orthosis

Provides proper support and immobilization of the wrist. Allows unobstructed motion of the thumb and fingers.

Suggested L-Code: L3984

### Indications

- Post cast fracture immobilization
- Wrist fractures and sprains
- Cumulative Trauma Disorder
- Carpal Tunnel Syndrome
- Post tendon transfer support
- Postoperative support after removal of internal or external fixation devices



## Wrist Hand Thumb Orthosis

Provides proper support and immobilization of the wrist and thumb. Allows unobstructed motion of the fingers.

Suggested L-Code: L3984

### Indications

- Post cast fracture immobilization
- Wrist fractures and sprains
- Cumulative Trauma Disorder
- Carpal Tunnel Syndrome
- Post tendon transfer support
- Postoperative support after removal of internal or external fixation devices
- Gamekeeper's thumb
- DeQuervain's Syndrome



Size	M-P Diameter*	Wrist Circum.	Forearm Circum.	Length of Brace	Wrist Hand product number	Wrist Hand Thumb product number
Small	2¾ - 3¼	6 - 7	8 - 10	7½	WH100 (L or R)	WHT100 (L or R)
Medium	3¼ - 3¾	7 - 8	10 - 12	7½	WH200 (L or R)	WHT200 (L or R)
Large	3¾ - 4	8 - 9	12 - 14	8	WH300 (L or R)	WHT300 (L or R)

**Southeast Medical Products, Inc**  
5524 Commerce Dr  
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)  
800.330.0890 (ph) • 800.851.9012 (fax)