

A grayscale X-ray image of a human knee joint, showing the femur (thigh bone) at the top, the tibia (shin bone) at the bottom, and the patella (kneecap) in the center. The joint space and surrounding structures are clearly visible.

**Southeast
Medical
Products, Inc**

***Specializing in
Orthopedics***

Since 1977

2011 Product Catalog

For over 30 years, Southeast Medical Products has provided the professional medical community with a great selection of quality products at fair and competitive prices.

Southeast Medical began in 1977 with a single product line, and today we represent dozens of manufacturers and warehouse thousands of inventory items, giving us the ability to offer a wide range of products and ship most orders the same day.

Southeast Medical Products prides itself on providing great service, too. We believe every customer and every order is important and we always strive to do whatever it takes to insure your satisfaction.

WE SINCERELY APPRECIATE YOUR BUSINESS!

BRANDS & MANUFACTURERS



...and many more.

Southeast Medical Products, Inc

5524 Commerce Dr

Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)

800.330.0890 (ph) • 800.851.9012 (fax)

ABDUCTION PILLOWS

Hip Abduction Pillow

All-foam construction, with concave sides to hold legs in position. Fire-retardant polyester foam with contact closure straps. Universal designs accommodate wide variation in patient sizes.

#79-90180 has 4" tear-off block on both ends of pillow for adjustment.

#79-90190 has removable, flannel lined, washable cover.



Part Number:	Width	Length	Size
79-90180	Bottom: 6-1/2" Top: 3"	27"	U
79-90173	Bottom: 12" Top: 6"	18"	S
79-90175	Bottom: 15" Top: 6"	23"	M
79-90177	Bottom: 18" Top: 6"	26"	L
79-90190	Bottom: 18" Top: 6"	18"	U

ACROMIOCLAVICULAR SUPPORT

Sturdy cotton/duck construction and felt shoulder straps provide additional comfort on shoulder area. Splint helps stabilize shoulder, humerus, and elbow following injuries. One size fits most patients.



Part Number:	Circ.	Size
79-85010	21"-42"	U

Measurement Information:

Chest circumference measurement taken at bottom of sternum.

AFO / DROP-FOOT SPLINTS

Pero-Max AFO

- The PEROMAX AFO limits plantar movement and provides a high energy return
- Constructed from a lightweight Carbon Fiber and contains an adjustable calf-bar, with two-calf pads, for larger sizes
- The supportive weight-bearing foot plate is rated for a maximum of 275 lbs.
- Flexible at the heel and forefoot
- The PEROMAX AFO provides excellent and dynamic gait



Indications:

Peroneal paralysis
Ankle-joint arthritis
Achilles post-operative treatment

<i>Height (cm):</i>	<i>Side:</i>	<i>Size:</i>	<i>Shoe Size:</i>	<i>Item#:</i>
38 cm	Left	Small	36-39	197P4/L1
38 cm	Right	Small	36-39	197P4/R1
39 cm	Left	Medium	39-42	197P4/L2
39 cm	Right	Medium	39-42	197P4/R2
40 cm	Left	Large	42-45	197P4/L3
40 cm	Right	Large	42-45	197P4/R3
41 cm	Left	X-Large	45-48	197P4/L4
41 cm	Right	X-Large	45-48	197P4/R4

Mfg suggested
HCPCS code; L-1951

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

AFO / DROP-FOOT SPLINTS

Matrix™

- Carbon composite
- Anterior shell is height-adjustable and heat moldable for customization and comfort
- Elongated strut is made of reinforced carbon composite for a more dynamic gait
- Anterior shell is padded reducing pressure on the tibial crest
- Custom-fit footplate
- Includes two washable pads for tibial crest relief



Mfg suggested HCPCS code; L-1932

Matrix MAX

- A rigid version of Matrix for more active and involved patients
- Sturdy strut is designed to offer stability to patients with weak quadriceps muscles
- Made of carbon and glass fibre composite
- Height-adjustable anterior shell allows optimal customization and comfort
- Includes two washable pads for the anterior shell
- Provides dynamic assist and reduced shearing for those with forefoot amputations



Mfg suggested HCPCS code; L-1932

Matrix Curve

- Optimal energy efficient gait
- Off-the-Shelf AFO with custom-fit calf band and footplate
- Lateral strut
- Open calcaneal design
- Slim, lightweight carbon design



Mfg suggested HCPCS code; L-1951

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

AFO / DROP-FOOT SPLINTS

Matrix™

Left	Right	Size	Height	Foot
FH200L-02	FH200R-02	S	15.5-16.5"/39.5-42cm	up to 25cm
FH200L-03	FH200R-03	M	16.0-17.0"/41.0-43cm	24-27cm
FH200L-04	FH200R-04	L	17.0-18.0"/43.0-46cm	25-29cm
FH200L-05	FH200R-05	XL	17.0-18.0"/43-46cm	26-31cm



Mfg suggested HCPCS code; L-1932

Matrix MAX

Left	Right	Size	Height	Foot
FH250L-02	FH250R-02	S	15.5-16.5"/39.5-42cm	up to 25cm
FH250L-03	FH250R-03	M	16.0-17.0"/41.0-43cm	24-27cm
FH250L-04	FH250R-04	L	17.0-18.0"/43.0-46cm	25-29cm
FH250L-05	FH250R-05	XL	17.0-18.0"/43-46cm	26-31cm



Mfg suggested HCPCS code; L-1932

Matrix Curve

Left	Right	Size	Height	Foot
FH275L-02	FH275R-02	S	15.5-16.5"/39.5-42cm	up to 25cm
FH275L-03	FH275R-03	M	16.0-17.0"/41.0-43cm	24-27cm
FH275L-04	FH275R-04	L	17.0-18.0"/43.0-46cm	25-29cm
FH275L-05	FH275R-05	XL	17.0-18.0"/43-46cm	26-31cm



Mfg suggested HCPCS code; L-1951

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

ANKLE SUPPORTS

Aircast Air-Stirrup® Ankle Brace

Since 1978, the Aircast Air-Stirrup Ankle Brace has been the "standard of care" for the functional management of ankle injuries and has been cited in over 100 medical journals for its superior performance in helping to heal ankle injuries. Each Aircast Ankle Brace features anatomically designed shells lined with the patented Duplex™ aircell system. This exclusive system incorporates two pre-inflated overlapping aircells, distal and proximal, that provide support and produce graduated compression during ambulation. The compression promotes efficient edema reduction in addition to helping accelerate rehabilitation. All Air-Stirrup Ankle Braces (except for the pediatric model) come with a patient guide providing information on brace application and optional rehabilitation exercises. The Standard (Large), Training (Medium), and Small size Ankle Brace is available individually or as part of an Ankle Sprain Care Kit™. The Pediatric Ankle Brace is designed for children 2 to 6 years of age and is ideal for ankle stabilization in hemophilia patients. Extra long straps are available "" contact Customer Service for more information.



Product Features and Benefits

- Semi-rigid, anatomically designed shells for protection, comfort, and the prevention of inversion/eversion
- Patented Duplex™ aircell system to enhance circulation and reduce swelling
- Pre-Inflated aircells for easy application
- Streamlined fit in shoes for early protected weight-bearing

Indication

- Acute injury
- Chronic instability

Left	Right	Size	Brace length
02AL	02AR	Standard (large)	10.5" (27cm)
02BL	02BR	Training (medium)	9.0" (23cm)
02CL	02CR	Small	8.7" (22cm)
02JL	02JR	Pediatric	6.0" (15cm)

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

ANKLE SUPPORTS

ASO Speed Lacer

Offers good ankle protection with ease of application.

Features & Benefits:

- **Speed lace closure:** allows easier and faster application
- **Stabilizing straps:** form figure-eight to protect and support ankle
- **Elastic cuff closure:** enhances support from stabilizing straps
- **Low profile:** will fit in any type of shoe
- **Bilateral design:** each size fits left or right foot
- **Breathable nylon boot:** enhanced breathability for long-term wear
- **US manufacturing:** higher quality control standards



Item Number	Size	Fits
223611	X-Small	Fits: 10" - 11"
223612	Small	Fits: 11" - 12"
223614	Medium	Fits: 12" - 13"
223615	Large	Fits: 13" - 14"
223616	X-Large	Fits: 14" - 15"
223617	XX-Large	Fits: 15" - 16"

Measure: Ankle Circumference around the heel



Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

ANKLE SUPPORTS

Bio Skin® TriLok™ Ankle Brace Prevent ankle sprains with the TriLok Ankle Control System

Whether you're looking for a brace to protect an ankle you've already injured or you're hoping to prevent an ankle injury, Bio Skin's® TriLok™ is the most effective brace for controlling inversion, eversion, and plantar flexion inversion. **No other brace offers more stability, comfort, and safety.**

- Inversion Control
- Plantar Flexion/Inversion Control
- Eversion Control
- High ankle sprains



Indications: Ankle sprains, chronic ankle instability, plantar fasciitis, posterior tibial tendon syndrome

Fit by shoe size	Small #53601	Medium #53602	Large #53603	Xlarge #53604
Mens	5-6.5	7-9.5	10-12.5	13-16
Womens	6-7.5	8-11	-	-

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

ANKLE SUPPORTS

Stabilized Ankle Support

Low profile design constructed with non-stretch ballistic nylon for comfort, performance and durability. Lace-up design and trimmable figure 8 strapping locks the calcaneus into place & controls inversion/eversion. A full circumferential elastic strap secures all strapping to provide superior tibia/fibula compression and stability. Removable medial/lateral stays for enhanced proprioception and stability. Fits easily into athletic or street shoe. Universal to fit right and left.



Part Number:	Circumference	Size
79-81352	10"- 11"	XS
79-81353	11"- 12"	S
79-81355	12"- 13"	M
79-81357	13"- 14"	L
79-81358	14"- 15"	XL
79-81359	15" - 16"	2XL
79-81359-10	16" - 17"	3XL

Measurement Information:

With athletic sock on, measure circumference of ankle from base of heel to around ankle.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

ANKLE SUPPORTS

Velocity Ankle Brace

The **Velocity** represents an evolution in DonJoy ankle brace technology, utilizing several new, proprietary design features, resulting in a new level of ankle brace performance and providing the foundation for successful ankle injury prevention and treatment. The **Velocity** is a low profile, lightweight ankle brace that prevents abnormal ankle inversion, eversion, and rotation while allowing natural, unrestricted dorsi- and plantar-flexion. The combination of soft goods with a rigid, hinged foot plate and calf cuff provides unsurpassed levels of control, fit, and support. Three **Velocity** versions are available, affording patients and athletes multiple ankle support levels.



Part Number:

81-1497x	Calf Cuff: EX - Standard, Left	EX=Extra Support
81-1498x	Calf Cuff: EX - Standard, Right	
81-1499x	Calf Cuff: EX - Wide, Left	
81-1500x	Calf Cuff: EX - Wide, Right	
11-1493-x-06000	Calf Cuff: MS - Standard, Left	MS=Moderate Support
11-1494-x-06000	Calf Cuff: MS - Standard, Right	
11-1495-x-06000	Calf Cuff: MS - Wide, Left	
11-1496-x-06000	Calf Cuff: MS - Wide, Right	
81-1489x	Calf Cuff: LS - Standard, Left	LS=Light Support
81-1490x	Calf Cuff: LS - Standard, Right	
11-1491-x-06000	Calf Cuff: LS - Wide, Left	
11-1492-x-06000	Calf Cuff: LS - Wide, Right	

Sizing Information: Foot Plate Sizes Available

x=2=Small	Men's Shoe: 6-8, Women's Shoe: 8-9.5
x=3=Medium	Men's Shoe: 8-12, Women's Shoe: 9.5-13.5
x=4=Large	Men's Shoe: 12 & up, Women's Shoe: 13.5 & up

Measuring Information:

Calf Cuff - measure calf circumference 8" from ground
Standard Calf: 7.25"-10" (18cm-25 cm)
Wide Calf: 10"-12.5" (25cm-32 cm)


Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

BOXER'S FRACTURE SPLINTS

TKO™ (The Knuckle Orthosis)



- Indicated for ligament injuries, Boxer's Fracture and other fractures of the metacarpals and proximal phalanges
- Straps are infinitely adjustable for fit and comfort (pain & swelling)
- Distal Buddy Strapping eliminates rotation of the finger
- Padded stays are malleable
- Reversible - can be used as a radial gutter splint



ComfortFORM™ Boxer's Splint



Durable lightweight fabric is **Lycra®** lined for breathability and patient comfort. The **ComfortFORM** Boxer's Splint is designed to eliminate the need for ulnar gutter casting or splinting. The malleable ulnar stay creates a dorsal block that maintains the 4th & 5th phalange at a 70° to 90° angle to keep fingers in full flexion following injury. The preformed palmer stay helps provide correct fit and proper support. Indications include Boxer's Fractures, flexor tendon and/or digital nerve damage to the 4th & 5th metacarpal phalange.



Boxer Splint

Anatomically designed to provide comfortable support, protection and immobilization of the 4th and 5th metacarpals

Indicated for treatment of acute or sub-acute fractures of the 4th and 5th metacarpals commonly referred to as "boxer fractures"

Constructed of lightweight low-density polyethylene shell with soft foam liner and palmer pad and contact closure straps



Extended Length

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

CANDO® EXERCISE BANDS

Exercise band therapy is the most commonly prescribed exercise modality for rehabilitation, conditioning and training. Simply cut an appropriate length of band and begin your exercise. The exercise band can be used with any band or tubing accessory.



Low Powder Latex Bands		
	<u>50 yd roll, 1/box</u>	<u>6 yd roll, 1/box</u>
	<u>order #</u>	<u>order #</u>
Tan	10-5220	10-5210
Yellow	10-5221	10-5211
Red	10-5222	10-5212
Green	10-5223	10-5213
Blue	10-5224	10-5214
Black	10-5225	10-5215
Silver	10-5226	10-5216
Gold	10-5227	10-5217
Latex Free Bands		
	<u>50 yd roll, 1/box</u>	<u>6 yd roll, 1/box</u>
	<u>order #</u>	<u>order #</u>
Tan	10-5620	10-5610
Yellow	10-5621	10-5611
Red	10-5622	10-5612
Green	10-5623	10-5613
Blue	10-5624	10-5614
Black	10-5625	10-5615
Silver	10-5626	10-5616
Gold	10-5627	10-5617

Clinic and Hospital Packs - Bulk rolls are packaged in 6 yd (5.5 m) and 50 yd (45.7 m) dispenser boxes. The boxes are color-coded to resistance level and design-coded to indicate whether the box contains low-powder or no-latex Cando exercise band, also helpful for inventory control and storage.



Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

CASTING MATERIALS

Whatever brand you prefer, you can count on ***Southeast Medical Products*** to deliver great pricing, quick delivery & unparalleled service.



3M™ Scotchcast™ Plus Casting Tape

- Water-activated resin lets you control working time according to application needs.
- High and early strength allows for weight bearing in 20 minutes and provides a durable cast with less chance of breakdown and the need for costly repairs or replacements.
- Lightweight for easy mobility.
- Stretches to conform well to contours for a more comfortable, better fitting cast.
- Porous weave allows the skin to breathe.
- Available in 11 colors to maximize patient satisfaction.



Delta-Lite® Conformable Fiberglass Cast Tape



- ✓ Exceptional conformability
- ✓ Three-dimensional stretch substrate
- ✓ Superior strength
- ✓ Weight bearing in approximately 20 minutes
- ✓ Excellent end lay down
- ✓ Reduces tucks or folds
- ✓ Superior lamination
- ✓ Easy open pouch
- ✓ 5-7 minute set time
- ✓ Available in ten popular colors

**Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839**

**407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)**

CAST REMOVAL SYSTEMS

American Orthopedic™

Power & performance at the right price



A/O Cast Saw
#0295-200

Specially designed for rugged use; perfect tool with power and endurance. Balance of weight and power helps reduce user fatigue. Internally lubricated for longer service. Switch is designed and placed so that it will not break if dropped. Hex drive provides 6 blade rotations extending blade life by 50% compared to similar competitive models. Length 12.5" Weight 3 lbs 13oz.

The A/O Cast Saw is compatible with your choice of vacuum



A/O Mobile Vacuum
#0295-400

AMERICAN ORTHOPAEDIC™ Vacuum's quiet operation eases patient fears through reduction of noise. However, it is powerful enough to pick up even fine particles. M-PACT vacuums are the only vacuums that incorporate a unique filter cartridge thus creating a safer environment. Once the dust enters the filter it is trapped. Other vacuums use paper bags that are messy to clean. The remote on/off switch gives the ability to use the vacuum alone or with the saw.



A/O Portable Vacuum
#0295-300

AMERICAN ORTHOPAEDIC™ Portable Vacuum's unique filtration system collects cast dust in a disposable filter cartridge- replacement is quick, easy and clean with no more messy bags. The convenient storage within the self-contained unit has space for cast spreader, scissors, replacement blades and the filter storage caps. The lightweight construction provides easy mobility. Furthermore, the portable vacuum is versatile as it attaches easily to most cast cutters.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

CAST REMOVAL SYSTEMS

Turbo Care

The Gold Standard in Cast Removal Systems



TurboCare Cast Saw w/Power supply
#TCCSS1

TURBOCARE® Cast Saw. Designed to reduce noise to a minimum results in less emotional trauma for patients and less hearing trauma for the clinician. The innovative motor plus saw design delivers high power in a smaller, lighter package. In addition, the unique design enhances the saw's balance and ergonomics for more user comfort. Plus automated multi-speed controls with "turbo-mode" for increased power.




TurboCare Vacuum
#TCCVS1

TURBOCARE® VACUUM. The vacuum's quiet operation eases patient fears through reduction of noise. However, it is powerful enough to pick up even fine particles. M-PACT vacuums are the only vacuums that incorporate a unique filter cartridge thus creating a safer environment. Once the dust enters the filter it is trapped. Other vacuums use paper bags that are messy to clean. The remote on/off switch gives the ability to use the vacuum alone or with the saw.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

CAST WEDGES



CastWedges aid in realignment and accommodate swelling. When used after an initial injury or surgery, the CastWedge permits uni- or bi-valving the cast, preventing extreme discomfort due to swelling. The CastWedge will maintain a precise opening in the cast. For realignment, CastWedge allows for exact adjustments to casted fractures using the opening wedge cast adjusting technique. The opening wedge cast/fracture technique avoids the inconvenience and cost of replacing the original cast. CastWedge is lightweight, inexpensive and can be used with plaster or fiberglass casts. It is made of high-impact, non-porous, radiolucent material that will not interfere with X-rays.



The CastWedge starter set consists of 100 wedges, 20 of each size (10mm, 15mm, 20mm, 25mm and 30mm). Replacement wedges are available in packages of 20 wedges.



Order #990-00

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

CERVICAL

Aspen Collars

Aspen Cervical Collars are designed to restrict cervical spine flexion, extension and rotation to promote patient recovery. Two-piece polyethylene design enhances patient comfort and is easily adjusted with Velcro™ tabs. Anterior and posterior flex tabs on the front and back panels provide improved rigidity and durability. The large tracheal aperture provides functional access for the airway management and other essential procedures. Removable, foam laminated pads wick away moisture to ensure proper patient hygiene and comfort while preventing skin irritation. X-Ray, CT and MRI lucent. Innovative height adjustment technology on VISTA and VISTA TX offers six (6) collars in one. Collar Set includes extra set of replacement pads.



Vista® Collar



Vista® TX Collar



Aspen® Collar

Order and sizing information on back

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

CERVICAL



Vista® Collar

Part Number:	Description	Circ.	Size
79-83380	Vista Collar set	13"-19"	Adult
79-83370	Vista Collar	13"-19"	Adult
79-83371	Vista Collar repl pads		Adult
79-83377	Vista Collar back		Large



Vista® TX Collar

Part Number:	Description	Circ.	Size
79-83360	Vista TX Collar Set	13"-19"	Adult
79-83390	Vista TX Collar	13"-19"	Adult
79-83371	Vista Collar repl pads		Adult
79-83377	Vista Collar back		Large



Aspen® Collar

Part Number:	Description	Height/Age	Neck Circ.
79-83363	Adult Short-Collar set	2.25"	13"-22"
79-83365	Adult Regular- Collar set	3.0"	13"-22"
79-83367	Adult Tall- Collar set	3.75"	13"-22"
79-83368	Adult XTall- Collar set	4.5"	13"-22"
79-83283	Adult Short	2.25"	13"-22"
79-83285	Adult Regular	3"	13"-22"
79-83287	Adult Tall	3.75"	13"-22"
79-83288	Adult XTall	4.5"	13"-22"
79-83272	Infant Short	1-18 months	8"-12"
79-83273	Infant Regular	9-24 months	11"-13"
79-83275	Child Tot	1-3 years	11.5"-14"
79-83277	Child Short	2-5 years	12"-14.5"
79-83278	Child Regular	3-6 years	13"-16"
Accessories			
79-83295	Replacement Pads	Child	
79-83325	Replacement Pads	Adult	
79-83383	Back	S	
79-83387	Back	L	

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

Miami J® Collar

Patented and proven, the Miami J® meets the rehabilitation needs of adults. Based on Clinical Studies and the input of leading medical institutions, this ground-breaking design provides the best immobilization to optimize patient rehabilitation, protects skin integrity during extended wear, and offers comfort features that enhance compliance and help speed recovery. All of our collars are latex-free and compatible with all imaging modalities.



Miami Jr.®

Pediatric c-spine care is distinct from adult care. The reason is basic — children are built differently. The Miami Jr.® from Jerome is the breakthrough collar designed exclusively for children 12 and under. The Miami Jr. line spans the growth cycle, is sized according to the Broselow-Luten color-coding system and meets the developing anatomical requirements unique to children.

When you snap on the Miami JTO Thoracic Extension, you make a great collar even better. Designed in conjunction with leading spine surgeons, orthotists and nurses, the JTO is used exclusively with all Miami J collars when even greater stability at the cervical junctions is indicated.



**Miami JTO®
Thoracic Extension**

**Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839**

**407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)**

MIAMI COLLARS

PRODUCT INFORMATION

Effective March 1, 2005



Size	Miami J [®] Collar Sets (collar + replacement pads)	Miami J [®] Collar	Miami J [®] Replacement Pads
Specialty Sizes:			
Super Short**	MJR-200S	MJ-200S	MJP-100
XS	MJR-250	MJ-250	MJP-250
Stout	MJR-200L	MJ-200L	MJP-200L
Standard Sizes:			
Small	MJR-300	MJ-300	MJP-100
Medium	MJR-400	MJ-400	MJP-100
Large	MJR-500	MJ-500	MJP-100

Size	Miami J [®] Front Panels	Miami J [®] Back Panels
Specialty Sizes:		
Super Short**	MJ-200S FRT	MJ-400 BK
XS	MJ-250 FRT	MJ-250 BK
Stout	MJ-200L FRT	MJ-200L BK
Standard Sizes:		
Small	MJ-300 FRT	MJ-300 BK
Medium	MJ-400 FRT	MJ-400 BK
Large	MJ-500 FRT	MJ-500 BK

Miami J[®] Front and Back Panels may be interchanged to accommodate patient anatomy, e.g., a Small Front may be used with a Large Back.

**Super Short collars are packaged with Medium Back Panels



Size	Miami Jr. [®] Collar Sets (collar + replacement pads)	Miami Jr. [®] Collar	Miami Jr. [®] Replacement Pads
0 – 6 months*	MJR-P0	MJ-P0	MJP-P0
6 months – 2 years	MJR-P1	MJ-P1	MJP-P1
2 – 6 years	MJR-P2	MJ-P2	MJP-P2
6 – 12 years	MJR-P3	MJ-P3	MJP-P3

* P0 collars are packaged with P1 Back Panels



Miami Occian Collar Back	ACB-100
--------------------------	---------

Southeast Medical Products, Inc
 407.851.5848/800.330.0890 (ph)
 407.851.9012/800.851.9012 (fax)
 orders@southeastmedical.com

CERVICAL

Philadelphia Tracheotomy® Collar

Anterior and posterior rigid plastic reinforcement supports limit movement. Total cervical arch support maintains neutral alignment of the cervical spine.

- Water-resistant for bathing and aqua-therapy
 - Easily customized for individual requirements
 - Latex-free, non-toxic and hypoallergenic
- Plastazote® foam reduces skin irritation



	Neck Circ ↓	1.25"	1.75"	2.25"	3.25"	4.25"	5.25"
Infant	6-8"	PHP-TI					
Pedi	8-11"		PHP-TP				
Xsmall	Up to 10"			PHP-T2XS	PHP-T3XS	PHP-T4XS	PHP-T5XS
Small	10-13"			PHP-T2S	PHP-T3S	PHP-T4S	PHP-T5S
Medium	13-16"			PHP-T2M	PHP-T3M	PHP-T4M	PHP-T5M
Large	16-19"			PHP-T2L	PHP-T3L	PHP-T4L	PHP-T5L
Xlarge	19" +			PHP-T2XL	PHP-T3XL	PHP-T4XL	PHP-T5XL

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

C E R V I C A L

Pneu-trac® Cervical Traction Collar

A comfortable, lightweight and easy-to-use alternative to traditional over-door traction, the Pneu-trac® provides cervical stretch through an inflatable collar. Clients can vary air pressure appropriate to the amount of cervical distraction required by using the hand-held inflation bulb with an air release valve. Made of semi-rigid polyethylene. Comes in three sizes.



- Pneumatic distraction control
- Patient adjustable
- Lightweight, semi-rigid polyethylene construction
- Two piece construction allows for varying height adjustments
- Velcro closure

Re-order #	Size	Neck circ.
A26700	Small	12¼ - 14¼"
A26701	Medium	14½ - 17"
A26702	Large	17 - 20"

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

CERVICAL

EMT'S[®] CHOICE COLLAR

Multi-Height Adjustments fit most adults. Essentially flat, one-piece design enables efficient storage where space is limited.

FEATURES & BENEFITS

- Restricts flexion, extension and rotation to immobilize the cervical spine.
- Injection molded polyethylene plastic body for improved rigidity and durability.
- Latex-free, non-toxic and hypoallergenic foam reduces skin irritation.
- Oversized trachea opening for quick access to carotid pulse monitoring and emergency tracheotomies.
- Pre-molded level chin support maintains patient in neutral alignment.
- Mandible support adapts to individual patient's unique jaw line.
- Fluorescent locking clips that positively lock into position.
- Open rear ventilation panel for fluid drainage, airflow and cervical spine palpitation.
- Water-resistant for use in adverse conditions.
- X-ray, CT and MRI lucent.



U.S. Patent 5,622,529 / 5,795,315 / 6,071,255 / 6,090,058 / D393,718



Life Without Limitations

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

CERVICAL

EMT'S® CHOICE COLLAR

MEASUREMENT & SIZING

Proper sizing is important for patient immobilization and comfort. Always maintain the patient's head in neutral alignment.



Height:

Measure the vertical distance from the top of the patient's shoulder to the tip of the chin (Fig. A).

The measurement corresponds to the HEIGHT of the Collar.

If the patient's measurement falls between two consecutive sizes, apply the smaller size first.



Circumference:

Measure around the patient's neck (Fig. B).

The measurement corresponds to the CIRCUMFERENCE of the Collar.

APPLICATION

Proper application of the Collar is as important as proper sizing for immobilization and patient comfort. The one-piece design consists of a single piece which is packaged as a unit.



Size Collar:

First, size the patient using your fingers. Then, find the **SIZING LINE** on the Collar and position your lowest finger on the line. Select the size closest to the height of your fingers.



Lock Size:

Position the Collar with the selected sizes. Press the lock into the Collar until secured. Verify by attempting to readjust. If movement occurs, repeat procedure until the Collar is secure.



Adjust Size:

To unlock, simply push the lock to return the lock to its unsecured position. Determine the proper height for the patient and press the lock into the selected size. Verify by attempting to readjust. If movement occurs, repeat procedure until Collar is secure.



After selecting the properly sized Collar, apply the Collar to the patient's neck. The chin should be secured by the patented molded chin support.

If the patient is in a supine position, slide the **BACK** part of the Collar behind the patient's neck until the Velcro® strap is visible.



Before attaching the Velcro® strap, make certain the patient's head and the Collar are in neutral alignment. If the chin is not sufficiently supported, tighten Collar. If additional tightening results in patient extension, select smaller size.

SIZE CHART & PRODUCT CODES

Size
ADULT

Product Code
PHP-EMT

Height
UNIVERSAL

Circumference
11" – 23"



Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

Original Knee Strap

For two decades, orthopedic surgeons, rheumatologists, podiatrists, family-practice physicians, physical therapists and their patients, professional athletes, and active individuals have praised our Knee Strap. Why? Shortly after initial usage, it provides many individuals with a marked reduction in the pain and discomfort caused by various knee disorders -- chondromalacia patella, Osgood Schlatter's Disease, runner's and jumper's knee, some forms of arthritis, and knee degeneration -- while, at the same time, allows them to return to their normal activities with increased mobility.



Cho-Pat's Original Knee Strap stabilizes and tightens up on the kneecap mechanism by applying pressure upon the patellar tendon below the kneecap. This compression reduces or eliminates inflammation and helps prevent knees from giving out.

Re-order # Black / Tan	Size	Circumference below patella
CP-12 / CP-13	Extra Small	Less than 10 in.; Less than 25 cm.
CP-14 / CP-15	Small	10-12.5 in.; 25-32 cm.
CP-16 / CP-17	Medium	12.5-14.5 in.; 32-37 cm.
CP-18 / CP-19	Large	14.5-16.5 in.; 37-42 cm.
CP-20 / CP-21	Extra Large	16.5-18 in.; 42-46 cm.
CP-11 / CP-22	Super Large	18-21 in.; 46-53 cm.

Dual Action Knee Strap

Cho-Pat's newly patented Dual Action Knee Strap provides an extra dimension of relief for painful and weakened knees.

The strap applies pressure upon the patellar tendon below the kneecap to stabilize and tighten up the kneecap mechanism, which improves patellar tracking and elevation and reduces patellar subluxation. It adds further strengthening to the kneecap mechanism by applying pressure on the tendon above the kneecap as well. This tends to reduce the forces of the quadriceps on the patella tendon and erosion of the undersurface of the kneecap due to a possible misalignment of the quadriceps (Q angle).



The device's design and construction apply constant dynamic forces to the surrounding areas of the knee reducing the likelihood of overuse syndromes. At the same time, the strap allows for full mobility while providing maximum benefits.

Re-order#	Size	Circumference @ mid-knee
CP-01	Extra Small	10-12 in.; 25-30 cm
CP-02	Small	12-14 in.; 30-36 cm.
CP-03	Medium	14-16 in.; 36-41 cm.
CP-04	Large	16-18 in.; 41-46 cm
CP-05	Extra Large	18-20 in.; 46-51 cm.
CP-07	XX Large	20-22 in.; 51-56 cm.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

Counter Force Knee Wrap

Cho-Pat's Counter-Force Knee Wrap is a versatile and dynamic approach to the pain and discomfort associated with knee degeneration resulting from old injuries, arthritis, wear and tear, or normal aging.

Constructed of heavy-duty, stretchable, fabric-covered neoprene, the Counter-Force Knee Wrap provides support and warmth to the knee helping to alleviate patellar pain while still allowing the user full mobility. Two adjustable straps on the device with hook and loop fasteners make it possible to adjust the wrap for additional kneecap support.



Re-order#	Size	Circumference @ mid-knee
CP-91	Extra Small	12-13.5 in.; 30-34 cm.
CP-92	Small	13-14.5 in.; 33-37 cm.
CP-93	Medium	14-15.5 in.; 36-39 cm.
CP-94	Large	15-16.5 in.; 38-42 cm.
CP-95	Extra Large	16-17.5 in.; 41-44 cm.
CP-96	XX Large	17-18.5 in.; 43-47 cm.
CP-97	Super Large	18-19.5 in.; 46-50 cm.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

CUBITAL TUNNEL SPLINTS

Cubital Comfort™ Brace

- Ideal for treatment of Cubital Tunnel Syndrome
- Two aluminum malleable stays positioned inside the elbow prevent excessive elbow flexion
- Aluminum stays are pre-bent at 20 degrees of full flexion
- Two stockinette arm bands aid in arm positioning



Ordering information

Short < 11 1/2"	Long > 11 1/2"
#3847S	#3847L

Measure length from mid forearm to mid humerus

PDAC Approved L3762

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

CUBITAL TUNNEL SPLINTS



XR

- Innovative support that immobilizes and positions the elbow at 122 ° to help reduce pain and numbness associated with Cubital Tunnel Syndrome
- Indicated for Cubital Tunnel Syndrome, post-surgical application or post-cast removal
- Constructed of rigid A.B.S. material with medial contouring on shell and plush cloth liner (replacement liners also available)



Ventral

- Unique anterior design immobilizes and positions elbow at 122° while allowing easy access to dressing
- Indicated for Cubital Tunnel Syndrome or post-op application
- Constructed of high density polyethylene shell with plush cloth liner

	S/M	M/L	L/XL
XR	CF37-2191	CF37-2192	CF37-2193
Ventral	CF37-2291	CF37-2292	CF37-2293
	Most females & smaller framed males	Most males & larger framed females	Large framed males & XL framed females

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

CUBITAL TUNNEL SPLINTS

IMAK Elbow Wrap

Provides comfort and pain relief while restricting painful movements during rest and sleep. Removable splint allows maximum immobilization with the stay inserted, or gentle support with the stay removed. Two adjustable straps for easy change of size and compression. Constructed with breathable foam covered with soft poly flannel for support and warmth. Latex free. Ideal for relief of cubital tunnel syndrome, Medial epicondylitis (Golfer's elbow), Lateral epicondylitis (Tennis elbow), and other elbow conditions.



Ordering Information

Measurement taken at forearm circumference

Part Number	Measurement	Size
79-81190	≤ 11"	Universal
79-81210	≥ 11"	X-Large / Universal

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

D M E - CRUTCHES/CANES/WALKERS



Lightweight aluminum crutches feature a double extruded center tube that provides additional strength to weight-bearing area. Push-pin adjustment makes it easy to adjust crutch height in 1" increments. Non-skid, jumbo size vinyl contoured tips provide excellent traction.

Item #	Size	Patient height
10400	Adult	5'2" – 5'10"
10401	Youth	4'6" – 5'2"
10402	Tall	5'10" – 6'6"
10416	Pediatric	4' – 4'6"



Adjustable aluminum cane features ergonomically designed handle with soft foam grip for comfort and security. Height adjusts from 30" to 39". Manufactured with sturdy, 1" diameter anodized extruded aluminum tubing.

Item #	Color	Adjustable Height
10306	Black	30" – 39"
10307	Bronze	30" – 39"



Quad Cane

Item #	Color	Base
10300	Chrome	Large
10301	Chrome	Small



Aluminum walker features single-button release to allow quick, one-hand folding. Cross brace is located near hand grips to allow for full stride. Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight. Includes rear glide caps, (item # 10107) and glide covers (item #10107C) allowing use on all surfaces.

Item #	Size	Adjustable Height
10202	Adult	32" – 39"
10203	Youth	25" – 32"
10212	Adult w/wheels	32" – 39"



5" wheels, item#10109



Hemi-Walker, item #10240

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

D M E - CRUTCHES/CANE/WALKERS

K. D. (knock down) Universal Aluminum Crutches



Do you think this is a tall crutch and a youth crutch? Well, you're right. But it's not two separate crutches, it is the **same crutch!** Reduce your inventory needs and free up space in your stock room with K.D. Universal Aluminum Crutches available from Southeast Medical Products, Inc.

- Universal height adjustment for use by a child, adult and tall adult. (individuals 4'6" to 6'6")
- New Euro-style clips and push pin adjustment- make it easy to adjust hand grip.
- Lighter than wood, the aluminum crutch provides stability and durability.
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction.
- Comfortable, durable underarm pad and hand grip.
- Weight capacity: 300 lbs
- Boxed (8" x 4.5" x 21.5") for easy storage and to take up less space in stock room.
- Limited Lifetime Warranty
- **Order #RTL10433**



Prefer standard crutches? We have those too!

Item #	Size	Patient height
10416	Pediatric	4' - 4'6"
10401	Youth	4'6" - 5'2"
10400	Adult	5'2" - 5'10"
10402	Tall	5'10" - 6'6"

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

DME - CRUTCHES / HEAVY DUTY

ReBound Spring Assisted Crutch



Patient Type

- Adolescents or adults from 4'7" to 6'10"
- **Weight capacity of 400lbs per crutch (800lbs per set)**

Features and Benefits

- Shock absorber/Spring assist
 - Absorbs impact and returns captured (negative) energy, helping the user to move forward
 - Eliminates painful impact in underarm area
 - Helps reduce pain in back, shoulders, arms, and wrists
- Ergonomic handle
 - Keeps hand and wrist in neutral position
 - Allows proper blood flow and nerve conduction
 - Helps reduce carpal tunnel syndrome and wrist tendonitis
- Folding design
 - Convenient to fold, store, and travel with
- Underarm cradle
 - Extended area creates more comfortable feel
 - Helps reduce underarm soreness and irritation
 - Eliminates axillary artery stenosis and thrombosis

Item #	Size	Patient height
11-0252-1	Short	4'7" – 5'7"
11-0252-2	Tall	5'7" – 6'10"

Recommended HCPCS Code: E0114

Bariatric Steel Crutch

- **500 lb. Weight Capacity**
- Durable, steel construction.
- Push-buttons easily adjust height in 1" increments.
- Comfortable, contoured underarm design.
- Easy wing nut adjustments for hand grips.
- Large, non-skid jumbo size vinyl contoured tips.
- Limited Lifetime Warranty.



Item #	Size	Patient height
10406	Adult	5'2" – 5'10"
10408	Tall	5'10" – 6'6"

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

DME - FOREARM CRUTCHES



- Vinyl coated, ergonomically contoured arm cuffs molded for comfort and stability.
- Leg and forearm sections adjust independently for optimal sizing.
- Vinyl hand grips are comfortable and durable.
- Extra large tips for added stability.
- Ortho K handle is hard plastic.
- Limited Lifetime Warranty.

Item #	Description	Patient height	Handle height	Weight limit
10403	Adult	5' 0" – 6' 2"	28" – 37"	300 lbs
10404	Youth	4' 2" – 5' 2"	21.5" – 30.5"	300 lbs
10405	Tall	5' 10" – 6' 6"	33" – 42"	300 lbs
10407B	Child, Blue	3' 2" – 4' 6"	14" – 21"	250 lbs
10407R	Child, Red	3' 2" – 4' 6"	14" – 21"	250 lbs

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

DME - ROLLATORS

Aluminum Rollator w/Fold Up and Removable Back Support, Padded Seat, 8" Casters



- Large 8" casters are ideal for indoor and outdoor use.
- Removable, hinged, padded backrest can be folded up or down as necessary.
- Handles are adjustable in height.
- Easy to use loop locks.
- Comes with padded seat with zippered pouch under seat.
- Comes standard with pouch.
- Optional basket (item # 9501S72606).
- Limited Lifetime Warranty.
- Available in Red (shown), Blue or Black

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

DME - SPECIALTY WALKERS

DV8 Steerable Knee Walker model #790



- Knee walker can be steered for increased maneuverability.
- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot.
- Deluxe lever brake ensures safety.
- Tool free height adjustment.
- 8" casters are ideal for indoor/outdoor use.
- Limited Lifetime Warranty.
- Lifetime Warranty on brake cable.

Economy Knee Walker model #780



- Tool free change from left leg to right leg use.
- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot.
- Deluxe lever brake ensures safety.
- Tool free height adjustment.
- 8" casters are ideal for indoor/outdoor use.
- 300 lb. Weight Capacity on 780 and 780J.
- 400 lb. Weight Capacity on 780HD.
- Limited Lifetime Warranty.
- Lifetime Warranty on brake cable.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

DORSAL CARPAL TUNNEL

The RCAI Dorsal Carpal Tunnel Splint is made of a lightweight, durable and heat moldable thermoplastic for patient customization. The narrow palmar roll supports the hand while fingers remain free to grasp.



The Dorsal Carpal Tunnel Splint limits wrist flexion and extension during repetitive hand motion. It can also be used as a wrist extension splint to support a weak or flaccid wrist.

The breathable foam liner wicks moisture away from the skin.

Size	Measure width of palm	Model # Right	Model # Left
Ex Small	2" - 2 1/2" (5.1 - 6.4 cm)	44DCT-ES-R	44DCT-ES-L
Small	2 1/2" - 3" (6.4 - 7.6 cm)	44DCT-S-R	44DCT-S-L
Medium	3" - 3 1/2" (7.6 - 8.9 cm)	44DCT-M-R	44DCT-M-L
Large	3 1/2" - 4" (8.9 - 10.2 cm)	44DCT-L-R	44DCT-L-L
X-Large	4" - 4 1/2" (10.2 - 11.4 cm)	44DCT-XL-R	44DCT-XL-L

EDEMA GLOVES

Soft, comfortable compression gloves for edema patients.
Available in 3/4" or full finger length. Color: Champagne.



Part **Description**

Part	Description		
902XSL	7"	X-Small	Left
902XSR	7"	X-Small	Right
902SL	8"	Small	Left
902SR	8"	Small	Right
902ML	9"	Medium	Left
902MR	9"	Medium	Right
902LL	10"	Large	Left
902LR	10"	Large	Right

Part	Description		
903XSL	7"	X-Small	Left
903XSR	7"	X-Small	Right
903SL	8"	Small	Left
903SR	8"	Small	Right
903ML	9"	Medium	Left
903MR	9"	Medium	Right
903LL	10"	Large	Left
903LR	10"	Large	Right



Part	Description		
904SL	8"	Small	Left
904SR	8"	Small	Right
904ML	9"	Medium	Left
904MR	9"	Medium	Right
904LL	10"	Large	Left
904LR	10"	Large	Right

Specify RIGHT or LEFT. Measure around the hand, palm flat, around the knuckles excluding the thumb. Order larger size if between sizes.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

ELBOW BRACES, HINGED



ProCare ElbowRANGER® Motion Control Splint

For non-operative application, collateral ligament strains and ruptures, and ulnar nerve transposition. May be used following early cast removal to stabilize fractures, and for post-operative immobilization. Removable hand grip assists in the prevention of pronation/supination. Aluminum R.O.M. bilateral hinges allow adjustments in 30° increments with formed Kydex® shells. Universal application fits right or left arm.



DonJoy IROM Elbow

The IROM Elbow is ideal for fixed or controlled ROM in the treatment of stable supracondylar fractures, stable radial fractures, and lateral/medial epicondylitis.



Ossur Innovator X®

Benefiting from our proven, patented technologies, Innovator X, with extending struts, can help improve patient outcomes. With superior comfort, clean aesthetics and the most user-friendly hinge in post-op bracing, Innovator X was designed to increase patient satisfaction and compliance.



BREG T-Scope Elbow Brace

The T-Scope Elbow is ideal for fixed or controlled range of motion. This lightweight, low profile product provides telescoping forearm and humeral bars, range of motion control, and lock out ability at any angle. Additionally, a detachable shoulder sling is included.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

HEEL / ELBOW PRESSURE RELIEF

Heelbo® Heel and Elbow Protector

Helps reduce shear and friction to heels and elbows while allowing skin to breathe.

- Flexible stretch weave fits comfortably and slips on easily
- Color-coded system allows you to determine the size quickly
- Sew-in contoured foam pad eliminates rolling and binding
- Machine washable; protectors retain shape after repeated use
- Available in 5 color-coded sizes



Cradle-Lite™ Elbow and Heel Protector

- Helps eliminate shearing and constriction
- Comfor™ foam inner pad maximizes skin aeration and minimizes pressure
- No irritating straps or ties
- Available in four sizes to assure proper fit
- Attractive color-coded sizing
- Comfortable circular knit sleeve
- Latex Free
- Available in 4 color-coded sizes



Flotation Elbow Pad

Visco-elastic foam distributes pressure evenly.

Knitted cotton, Orlon® and Lycra® sleeve with a visco-elastic foam cushion. Latex free. Protects the elbow from pressure sores or to lessen nerve hypersensitivity. Snug-fitting sleeve reduces migration of the pad on the arm. Sold in pairs.

To size, measure circumference at elbow crease.

Available in two sizes; Regular (9 ½" - 10 ½") and Large (10 ½" - 12")
measure elbow circumference



Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

ELECTROTHERAPY

Intelect® Portable Electrotherapy

For symptomatic relief and management of chronic pain and acute post-surgical or post-trauma pain.



- Dual, independent channels
- Adjustable parameters
- 5 stimulation modes (Burst, Normal, MRW, SD, Bi-pulse)
- Patient compliance timer
- Digital LCD display

#77712 Intelect® TENS Digital



- Dual, independent channels
- Adjustable parameters
- 3 stimulation modes (Burst, Normal, Modulation)
- Flip cover to protect top controls

#77600 Intelect® TENS Standard



- Dual, independent channels
- Adjustable parameters
- 3 stimulation modes (Burst, Normal, Modulation)

#77710 Intelect® TENS Basic

DURA-STICK® II
SELF-ADHESIVE ELECTRODES



Our most popular and best selling electrode, DURA-STICK® II electrodes provide comfortable and effective stimulation through the unique silver printed conductor.

Flexible, yet sturdy construction offers increased durability and ease of removal. Soft and conformable lead wire connectors offer superior grip to the lead wire pin.

All devices come with lead wires, 9-volt battery, DURA-STICK® II Self-Adhesive Electrodes and a convenient carrying case.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)


FINGER SPLINTS



STAX SPLINTS are designed to support distal joint of finger in extension while permitting unrestricted movement of proximal interphalangeal joint. Available in starter kit and individual sizes (sold each).



Digital splints in 11 pre-cut sizes meet most splinting needs. Available in starter kit and individual sizes (sold by the dozen)



Ceramic finish malleable aluminum strips with radius ends that can be cut to length needed. Available in blue, soft, open cell foam and white, firm, non-absorbent, closed cell foam. Sold by the dozen.



SPECIALTY KIT. An assortment of 45 popular padded specialty splints in a plastic organizer.



COT Aluminum splint folds evenly over the end of finger to provide stabilization and protection.



PROTECTOR Aluminum splint provides end protection to injured fingers.



SPOON Aluminum splint curves upward at tip to protect end of finger.



FROG Aluminum splint folds to hold finger in proper alignment.



BASEBALL Aluminum splint folds over middle of finger and easily adjusts and molds to finger.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

FINGER SPLINTS



Finger Splints stabilize and align the PIP and DIP joints. With a simple turn of the splint, each Oval-8 Finger Splint can be used to treat tendon or ligament injuries, protect healing fractures, support unstable joints, prevent contractures, and more.

Treat **six** of the most common finger conditions with just one splint.



Correct lateral instability



Prevent PIP hyperextension



Protect mallet finger injuries




Immobilize with two splints for protection



Prevent formation of Boutonniere deformities



Limit motion to prevent trigger finger



A fit for every finger

Available in sizes 2 to 15, each Oval-8 fits both a whole and a half size depending on which end is applied first. Our Sizing Set makes it easy to choose the correct size. Splints can also be adjusted with a Precision Spot Heat Gun (#P7500) for a perfect fit every time.



Oval-8 Kit - #P1008-K

With 44 splints in sizes 2 to 15, (the majority in the most commonly used sizes), the Kit provides everything needed to size and dispense the splint in seconds. Refills are available in packages of 5. Order #P1008-xx-5



Sizing Set – #P1008-S

Includes one of each size, 2 to 15, on an easy-to-open ring. Determine the correct size and order the splints as needed or dispense the splint immediately from the Oval-8 Kit.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

FINGER SPLINTS



3pp™ ThumSling™

Gently hugs the thumb and applies light compression to reduce grinding and pain at the CMC (basal) joint of the thumb.



Features

- Wrap around sling repositions the CMC joint at the base of the thumb so you can perform daily activities with less pain or fear of dropping things.
- Easy to apply, this soft, flexible splint allows you to grip a golf club, hold a rake or grip a pen securely and comfortably.
- Lightweight and lined in soft, breathable foam, the ThumSling is not hot or bulky. Machine washable.
- Light Control

Sizing Information

Instructions: Lay hand flat and measure around the MP joints (large knuckles).

Item Number	Size	Fits
P2006	Small/Medium	6" to 7 1/4"
P2006	Medium/Large	7 1/4" to 9 1/4"

Choose right or left hand.

3pp™ Buddy Loop™

Eliminates messy, sticky tape and prickly straps! Buddy Loops are the number one choice of physicians and athletes to protect finger injuries at rest and during play.



Features

- Fast "loop and wrap" design makes it easy to hold fingers together securely.
- Soft, foam-lined material provides a light cushion between your fingers for comfort and control.
- Use one or two Buddy Loops as needed for control. For larger hands or longer fingers, we recommend a wider one-inch strap at the base of the fingers used in combination with a half-inch strap near the fingertips.
- Machine washable.
- Light Control

Sizing Information

Instructions: 5" long loops fit most fingers. Available in three widths - 1/2" wide, 3/4" wide and 1" wide. Sold in packages of 5, 25, 50 or 100.

Item Number	Size
P1005	1/2" Buddy Loop
P1007	1" Buddy Loop
P1003	3/4" Buddy Loop

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

FINGER SPLINTS

LMB Light Flexion Finger Spring, #601

- High quality Wire-Foam™ construction for a customized fit
- Produces flexion of the PIP or DIP joints
- Finger pad is free allowing for function
- Capable of bending supple joint to approximately 80°
- May be worn on adjacent fingers
- Tension adjustable by careful bending - lighter than model #503


INDICATIONS: Crush injuries; Swan neck deformity



LMB Light Extension Finger Spring, #602

- High quality Wire-Foam™ construction for a customized fit
- Produces extensions of the PIP or DIP joints
- Capable of fully extending supple joint
- May be worn on adjacent fingers
- Tension adjustable by careful bending
- Lighter tension than models #501, #502, #507 and #508
- Force: Full flexion ≈ 1.75 lbs.; Full extension ≈ 1.0 lbs.

INDICATIONS: Crush injuries; Boutonniere deformity; Burns




Size	Measurement	FLEXION	EXTENSION
X-Small	1¼"	601AA	602AA
Small	1½"	601A	602A
Medium	1¾"	601B	602B
Large	2¼"	601C	602C
X-Large	2¾"	601D	602D

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

FOOT PRODUCTS



Rolyan® External Fixator Foot Support

A dynamic dorsiflexion assist for tibial fractures. Helps prevent foot drop and protects heel. Compatible with external fixation systems: including Hoffman, Kroner and Fisher. Elastic straps give dorsiflexion support and plantar flexion exercise. Can be positioned for inversion or eversion. Base washable.

External Fixator Replacement

Includes replacement foam, two elastic straps with snap catch, and one foam strap pad.

Item#	Description
A4981	Rolyan® External Fixator Foot Support
A4983	Rolyan® External Fixator Foot Support External Fixator Replacement

FOOT PRODUCTS / SPENCO



PolySorb Cross Trainer Insoles #38-034-xx
For all sports or athletic activity

For aerobics, court sports, running, weight training, racquet sports and walking. Special cushioning in both the heel strike zone and the forefoot area for outstanding performance. Helps turn any athletic shoe into a cross training shoe.



Slip-In Insoles #40-212-xx

Long lasting heel-to-toe comfort

Years of proven use and high quality, synonymous with the name SPENCO®, our slip-in insoles provide long lasting heel-to-toe comfort. Washable, trim to fit with scissors.



3/4 Length Orthotic Arch Supports #43-158-xx

Spenco® Orthotic Arch Supports are specifically designed to help control the position of the feet. They are shaped to help stabilize and balance the foot and to provide the extra comfort and support needed for optimum performance. This unique cushioning system helps absorb shock, reduce friction and improve overall foot comfort. Also available in full length.



Size Chart


xx	00	01	02	03	04	05	06
Women's	3/4	5/6	7/8	9/10	11/12		
Mens's			6/7	8/9	10/11	12/13	14/15

All Spenco products are available from Southeast Medical Products.

Southeast Medical Products, Inc
 5524 Commerce Dr
 Orlando FL 32839
 407.851.5848 (ph) • 407.851.9012 (fax)
 800.330.0890 (ph) • 800.851.9012 (fax)

FRACTURE BRACING

Maramed Humerus Fracture Brace



Indications

- Nonoperative treatment of closed humeral shaft fractures
- Postoperative treatment of open humeral shaft fractures
- Protection of humeral defects
- After removal of internal or external fixation
- Post biopsy
- Defects resulting from tumor growth
- Adjunct to internal fixation such as unreamed intermedullary nails

Item Number	Size Left/Right	Circumference	Medial Length	Lateral Length
HFB-100L or R	Small	8 - 11	5¼	6¼
HFB-200L or R	Medium	10 - 13	5½	6¾
HFB-300L or R	Large	12 - 15	6	7¼
HFB-400L or R	X-Large	14 - 17	6¼	7¾
HFB-500L or R	XX-Large	18 - 23	7¼	8¾

Suggested Billing Code: L3980

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

FRACTURE BRACING

Marmed Over the Shoulder-Universal Humerus Fracture Brace



Indications

- Nonoperative treatment of closed humeral shaft fractures
- Postoperative treatment of open humeral shaft fractures
- Protection of humeral defects
- After removal of internal or external fixation
- Post biopsy
- Defects resulting from tumor growth
- Adjunct to internal fixation such as unreamed intermedullary nails

Item Number	Size Universal	Circumference	Medial Length	Lateral Length
HFB-OS-100	Small	8 - 12½	6	11
HFB-OS-200	Medium	10 - 13¼	6½	12½
HFB-OS-300	Large	12 - 15	7	13½
HFB-OS-400	X-Large	14 - 16	7¼	14

Suggested Billing Code: L3980

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

HAND SPLINTS

Galveston Metacarpal Brace



- **ADVANTAGES:** Utilizes the concept of three point fixation, a better degree of correction than with the ulnar gutter splint, fully adjustable pads to customize the splint to various hand sizes and fracture types, and an adjustable strap that allows the maintenance of swelling reduction. It is radiotransparent, and allows unrestricted motion of IP, MP, and wrist joints.
- **INDICATIONS:** Can be applied to fractures of the 2nd through the 5th metacarpal on either hand by simply adjusting the crossbar and pads.
- **PRECAUTIONS:** The two-toned dorsal pad has a soft layer in contact with the skin, and a stiffer backing pad. As long as the soft pad is not fully compressed, pressure under the pad is below that shown to cause necrosis (the soft pad is only meant to be an indicator). As with other non-invasive methods of treatment, this splint does not prevent shortening, The Galveston Metacarpal Brace, as with any other functional bracing system, requires patient compliance for optimal results.
- **APPLICATION CONSIDERATIONS:** The brace should be applied following fracture reduction to maintain alignment. The patient should be allowed to adjust the brace to be snug, but not uncomfortable.

Sizing

Measure the Diameter of the Palm.



Order #	Description	Inches	Centimeters
GMB101S	Small	2.25"-3"	6-8cm
GMB101M	Medium	3"-3.5"	8-9cm
GMB101L	Large	3.5"-4"	9-10cm
GMB101XL	Extra Large	Over 4"	Over 10cm

A spare set of pads is included with each brace.
Extra pads & bars are available at an additional cost.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)


HEEL CUPS & PADS



Silicone Heel Cups

Molded medical-grade soft silicone to help reduce shock and vibration loads to joints, to address plantar fasciitis, heel spurs and achilles tendonitis. Softer center in heel region to help reduce strain to area sensitive to pain. Sold in pairs.

Number:	Male	Female	Size
79-81102	Children/youth	Children/youth	XS
79-81103	5-1/2 - 9-1/2	5 - 9	S/M
79-81107	10+	9-1/2+	L/XL



Aetrex Heel Spur Cradles were doctor designed to comfort heel pain. The cradle absorbs shock and cushions painful heels for immediate relief. Anti-Shox Spur Heel Cradle provides relief from heel spurs and plantar fasciitis.

Number:	Size
5075S	SML
5075M	MED
5075L	LRG



Bone Spur Pad

To reduce pain in the heel. Outer ring absorbs shock while inner ring provides an area to gently accommodate injury by containing soft tissue around heel. Relieves pressure on bone spur or in sensitive area of heel. Can help relieve lower back fatigue, heel spurs, plantar fasciitis and general heel discomfort.

Number:	Male	Female	Size
71301	4 - 6	5 - 7	SML
71302	6 1/2 - 10 1/2	7 1/2 - 11 1/2	MED
71303	11 +	12 +	LRG



Tuli's® Classic Heel Cups

Protect your heel in biomechanically-designed cradles that provide relief from heel pain.

Number:	Size
T10010	REG
T10020	LRG



Tuli's® Heavy Duty Heel Cups Equip your heel with the ultimate shock absorption required for high impact work or play.

Number:	Size
T10211	SML
T10212	MED
T10214	LRG



TuliGEL™ Heel Cups

Tuli's classic support in a lighter, softer but resilient gel cup.

Number:	Size
T10215	REG
T10216	LRG

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

HEEL / ELBOW PRESSURE RELIEF

Heelbo® Heel and Elbow Protector

Helps reduce shear and friction to heels and elbows while allowing skin to breathe.

- Flexible stretch weave fits comfortably and slips on easily
- Color-coded system allows you to determine the size quickly
- Sew-in contoured foam pad eliminates rolling and binding
- Machine washable; protectors retain shape after repeated use
- Available in 5 color-coded sizes



Cradle-Lite™ Elbow and Heel Protector

- Helps eliminate shearing and constriction
- Comfor™ foam inner pad maximizes skin aeration and minimizes pressure
- No irritating straps or ties
- Available in four sizes to assure proper fit
- Attractive color-coded sizing
- Comfortable circular knit sleeve
- Latex Free
- Available in 4 color-coded sizes



Flotation Elbow Pad

Visco-elastic foam distributes pressure evenly.

Knitted cotton, Orlon® and Lycra® sleeve with a visco-elastic foam cushion. Latex free. Protects the elbow from pressure sores or to lessen nerve hypersensitivity. Snug-fitting sleeve reduces migration of the pad on the arm. Sold in pairs.

To size, measure circumference at elbow crease.

Available in two sizes; Regular (9 ½" - 10 ½") and Large (10 ½" - 12")
measure elbow circumference



Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

L / S SUPPORTS - Warm/Form style



79-8250x

Low profile 9" elastic support with double side pull compression straps. Neoprene posterior panel accommodates ComfortFORM™ Moldable Insert #79-89120. Foam lumbar pad is included with 79-8250x style. No stays. Ideal for providing support and compression for strains & sprains of the low back and abdominal area.



79-8925x

9" neoprene/elastic construction. Low-profile design. Contact closure. Accommodates COMFORT-FORM Moldable Insert #79-89120.



79-8924x

Foam/pile, with 9" panel elastic construction and double side-pull elastic compression straps. Two sizes. Accommodates ComfortFORM™ moldable insert #79-89120. Ideal for providing support and compression for strains, sprains of the low back and abdominal area.



79-8249x

9" all-neoprene construction with double side-pull compression straps for additional support. 2-ply posterior pocket accommodates ComfortFORM™ Moldable Insert #79-89120. Ideal for providing support and compression for strains & sprains of the low back and abdominal area.



The ComfortFORM™ Moldable Insert can be heated with a heat gun or at 160° in a hot water pan, or hydrocollator and then formed against patient's back for exact contour and improved support. To be used with a compatible lumbosacral support to provide rigid support to the lower back.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839


407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

L / S SUPPORTS

Saunders Posture S'port®



The **Posture S'port** is unique because it anchors at the waist and not only discourages stooped shoulders, but assists in proper spinal alignment. The Posture S'port can be worn directly on the skin or over a t-shirt.



Features:

- Criss-crossed suspender like straps gently pull the shoulders back into the neutral position to decrease back and shoulder pain
- Soft cotton under arm pads provide a comfortable fit that won't pinch or bind
- Can be worn directly on skin or over clothing
- Adjustable elastic waistband that anchors the entire support to the waist

Small	22" - 30"
Medium	30" - 42"
Large	42" - 54"

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

MASSAGE THERAPY

Adapta® MC-100 Portable Massage Chair



Portability is a prerequisite for any successful massage therapist. Available in four colors: Royal Blue (#92019), Wedgewood (#92022), Maroon (#92020) and Mushroom (#92021).

Lightweight - 16 lbs. (7 kg)

UL listed at 300 lbs. (135 kg) working weight capacity

Durable and affordable

Shown with OPTIONAL Deluxe Carry Case, #93001.

**Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)**

PEDIATRICS - LOWER EXTREMITY



P E D I A T R I C S



Miami Jr



Patented design uniquely maintains airway and c-spine alignment of children during the growth transition from birth to adolescence.

Part Number:	Age
MJR-P0	0 - 6 mos
MJR-P1	6 mo - 2 yrs
MJR-P2	2 - 6 yrs
MJR-P3	6 - 12 yrs

ComfortFORM™ Wrist



Durable lightweight foam laminate is Lycra® lined for breathability and patient comfort. Preformed aluminum stay and loop/lock closure helps provide anatomically correct fit and proper support. Ideal for sprains, strains and control of carpal tunnel syndrome symptoms.

Left	Right	Wrist Circ.	Length	Size
79-87291	79-87281	Up to 4.5"	7.0"	Ped/XXS
79-87292	79-87282	4.5"-5.5"	7.5"	X-Small
79-87293	79-87283	5.5"-6.5"	7.5"	Small

Also available in pediatric sizes;

- Soft cervical collars
- Clavicle supports
- Shoulder immobilizers
- Finger/Thumb splints
- Neoprene Elbow sleeves
- Wrist/Thumb splints
- Posture Back supports
- And anything else we can get for you! Just let us know what you need.

Hinged Elbow



The Pediatric Universal Arm Brace has bilateral, lightweight stainless steel hinges, which provide increased flexion / extension elbow control from 0° to 120°. The precision length adjustment hinge arm assures a personalized fit and the unique forearm design provides protection of the olecranon process. The UAB is a cost effective alternative to casting for stable fractures and is indicated for elbow surgery, tendon or muscle repair, strains or sprains.

Humeral Fracture Brace



The Over the Shoulder Humeral Fracture Brace is a universal design and offers more support to humeral diaphyseal fractures with its proximal shoulder cap. Manufactured of trimmable polyethylene and lined with closed-cell foam, the Humeral Fracture Brace is lightweight, durable, and cannot absorb moisture.



Arm Slings

Easy slide buckle adjustment on webbing strap.

Part Number:		Size
79-99132	4-1/2" x 10"	Pediatric XS
79-99133	6-1/2" x 12-1/2"	S

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

PEDIATRICS - LOWER EXTREMITY



P E D I A T R I C S



Walking Boots



DonJoy Maxtrax
X-Smallsole
length is 8.75".



Bledsoe Wee Walker
Sole length 7.5" (AL185007) and
6" (AL185003)

Knee ImmoBs



Pile/foam laminate with moveable lateral/medial stays for exact positioning and easy adjustment. Alternating web straps provide dual adjustability, using a D-ring buckle and a hook tab. After fitting, knee immobilizer can be removed by simply releasing hook tab. Fits 29" maximum thigh circumference.

Part Number:	Length	Size
79-80600	8"	U
79-96870	10"	U
79-96871	12"	U

Aircast Ankle



Left	Right	Size
02JL	02JR	Ped, 6" brace length



Crutches

Ord #	Size	Fits
10416	Ped	4' - 4'6"
10401	Youth	4'6" - 5'2"

Cast Shoes

Darco Slimline



Ord #	Sole length	Children's Shoe
SLO1NP	5.5"	5 - 7.5
SLO2NP	6.5"	8 - 11.5
SLO3NP	7.5"	12 - 15

Also available in pediatric sizes;

- Neoprene Knee Sleeves
- Patella Stabilizers
- Ankle supports
- Knee Bracing
- Hip Bracing
- And anything else we can get for you! Just let us know what you need.

Rocker Bottom



Ord #	Sole length
0814-7298	4.5"
0814-7299	5.5"
0814-7300	6 3/8"

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

PEDIATRICS - HIP DYSPLASIA

KICKER" PAVLIK HARNESS

Biomechanically correct strap position for maximum effectiveness in treating hip dysplasia. We specialize in very small sizes suitable for the premie with hip dysplasia. Machine washable and dryable. Soft material in dark blue or light pink with self gripping fasteners.



Sizes:	
Premie	4-6 lbs.
Ex-Small	0-2 months
Small	0-3 months
Medium	3-6 months
Large	5-9 months

"CRUISER" HIP ABDUCTION BRACE

Polypropylene hip abduction brace for treating children from infants to 3 years of age with hip dysplasia. Holds hips in the correct abduction position, but flexible enough to allow walking. Can be used full time or for night wear only. Commonly used as a follow-up to the "Kicker" Pavlik harness treatment, as the child begins to walk.



Sizes:	Waist:
Ex-Small	14" - 16"
Small	16" - 18"
Medium	18" - 20"
Large	20" - 22"
Ex-Large	22" - 26"

Pediatric - Hip - Hip & Knee Abductor

The RCAI Pediatric Hip & Knee Abductor is constructed of Kydex® thigh cuffs and aluminum adjustment bar which controls scissoring of tight thigh adductors. Optional length turnbuckles maintain hip abduction between 20° - 60°. It offers a static position for bed or wheelchair or a reciprocating function for gait training / ambulation.



The liner is a washable synthetic sheepskin material.


Size	Measure Mid Thigh Circ.	Cuff Height	Model #
Pediatric	10 1/2 - 14" (26.7-35.6 cm)	4 1/2" (11.4cm)	16HR-P
Ex. Small	11 1/2 - 15" (29.2-38.1 cm)	5" (12.7cm)	16HR-ES
Small	13 1/2 - 16 1/2" (34.3-41.9 cm)	5" (12.7cm)	16HR-S

Suggested Billing Code L1650

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

PLANTAR FASCIITIS SPLINTS



procare™

Dorsiwedge Night Splint - Soft padded, contoured posterior shell provides passive dorsi-flexion to address plantar fasciitis and Achilles tendonitis. Fits right or left foot. Shipped with toe wedge for additional stretch on plantar fascia.



DARCO

The **Darco Night Splint** is designed to alleviate heel pain associated with plantar fasciitis and achilles tendonitis.

Special Features:

- Terry cloth lining is all natural, comfortable and moisture absorbent
- Non-slip tread on the sole
- Pneumatic dorsi-flexion air bladder provides infinite angle adjustment
- Maintains proper flexation when worn properly



BIRD & CRONIN INC.



Baker™ Plantar Fasciitis Night Splint

- *Designed to help minimize the pain and discomfort of plantar fasciitis*
- *Allows the patient to apply a consistent 90° angle stretch to the plantar fascia by holding the foot in gentle dorsiflexion*
- *Soft padded cover provides greater patient comfort and protection against irritation to the skin and breakdown*
- *Removable foam wedge insert is included for additional 5" stretch on the plantar fascia*
- *Fits either left or right foot*
- *Latex Free*

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

PLANTAR FASCIITIS SPLINTS



P.F.S.™ Plantar Fasciitis Stretch Splint

- *Designed to provide relief from pain and discomfort of plantar fasciitis and Achilles tendonitis*
- *Allows the patient to apply consistent stretch to the plantar fascia by holding the foot in gentle dorsiflexion*
- *Lightweight design for patient comfort while sleeping*
- *Slip resistant sole for greater patient safety*
- *Fits either left or right foot*
- *Latex Free*



Formfit Night Splint

- Adjustable flexion straps provide a gentle stretching of the plantar fascia and achilles tendon. Adjusts from 90 to 10 dorsi-flexion
- Patented 3-dimensional molding creates superior comfort and fit with padding where you need it
- Soft, lycra lined foam interior is plush, washable and breathable to eliminated heat build-up
- Lycra lined foam padded straps maintain comfort across the flexor tendons
- User-friendly center-release buckles simplify application
- Lightweight, low profile shell is sturdy and breathable



Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

PLANTAR FASCIITIS SPLINTS



Body Armor® Night Splint

The DARCO Body Armor® Night Splint is the only night splint on the market today that actively engages the windlass mechanism of the foot resulting in a specific and sustained stretch to the plantar fascia and a sustained low load stretch to the flexor tendons, Achilles tendon, and calf muscles.



- **Toe Strap** dorsiflexes the 1st MPJ to provide a specific and sustained stretch to the Plantar Fascia and a sustained low load stretch to the flexor tendons, Achilles tendon and calf muscles.
- **Loading the Windlass Mechanism** naturally brings the foot/ankle to a near 90 degree angle thus eliminating the “lever effect” and the resulting pressure on the apex of the ankle.
- **Elastic Connector Ring** allows for natural flexion of the foot and toes during sleep, returning to the dorsiflexed position when muscles relax.
- **Low Profile Design** allows for maximum user comfort and compliance.
- **One-size-fits-all** Design reduces inventory and ensures the right size is always on hand.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

PLANTAR FASCIITIS SPLINTS



AirForm Dorsal Night Splint

Performance Benefits:

- The Airform Night Splint provides an effective treatment for plantar fasciitis, achilles tendonitis, drop foot and post-static pain
- Its lightweight, low profile design makes it more comfortable. Better comfort means better compliance
- Slips on and adjusts easily
- Fits left or right foot
- Comes in 2 sizes

Indications:

- Plantar Fasciitis
- Achilles Tendonitis
- Drop Foot
- Post-Static Pain



The Strassburg Sock

For treatment of Plantar Fasciitis, heel spurs, Shin Splints and Achillies tendonitis.

Features:



- It is very easy to use.
- It does not lock the ankle joint into any one position thus allowing some movement between the various joints preventing possible joint compression pain.
- It can be safely worn for several hours without adverse effect when applied properly.
- It can be worn to bed, provides no significant restrictions and will not interrupt normal sleep if worn as directed.
- It is an effective device providing low cost relief of heel spur symptoms.
- It can reduce/eliminate the pain experienced with the first morning steps.
- It can, over a period of time, strengthen the arch of the foot.
- It can supplement and not interfere with the various traditional doctor's treatment.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

3x10=0 Full Extension Knee Device



The 3x10=0 Full Extension Knee Device is for the patient who lacks passive extension of the knee following ACL reconstruction, total knee replacement, knee surgery, and/or trauma. Upon first observation of knee extension lag, this prophylactic appliance allows patients to regain full extension through a safe and effective rehabilitation routine.

With the thin end under the thigh, the ankle is supported by the riser at the distal end of the foam cushion. The straps are drawn around the leg and the foam cushion. The straps are then tightened gradually to move the knee toward full extension.

The 3x10=0 Full Extension Knee Device is easier to use than the alternative day and night devices, therefore patient compliance is high. It can also be advantageous in overcoming long standing extension problems by simply lengthening the application time.

RESOLVE® HALO

Ossur's ReSolve Halo System is the industry's only MR-Safe halo system. ReSolve combines a glass composite ring with Ossur's proprietary ceramic-tipped skull pins.

It is the only halo clinically proven to eliminate the risk of dangerous scalp burns resulting from halo pin heating phenomena.

ReSolve isolates the patient from potential current induction in the [MRI](#). Not only does the Resolve system ensure MR-safety but clearer and more high-resolution brain images are now possible.



Features and Benefits

- The industry's only MR-Safe Halo System
- The only halo clinically proven to eliminate the risk of dangerous scalp burns resulting from halo pin heating phenomena
- Glass composite ring and ceramic-tipped pins are non-conductive and induce no dangerous current flow
- MR images are clear and artifact-free

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

RESOLVE® HALO

Part Number	Description
Resolve Nonconductive Halo Rings	
505300C	Resolve Closed Back Halo • 17-20" (43-51cm), Small
505400C	Resolve Closed Back Halo • 20-23" (51-59cm), Standard
505500C	Resolve Closed Back Halo • 23-26" (59-66cm), Large
510300C	Resolve Closed Back Halo, Ring Only • 17-20" (43-51cm), Small
510400C	Resolve Closed Back Halo, Ring Only • 20-23" (51-59cm), Standard
510500C	Resolve Closed Back Halo, Ring Only • 23-26" (59-66cm), Large
505300D	Resolve Open Back Halo • 17-20" (43-51cm), Small
505400D	Resolve Open Back Halo • 20-23" (51-59cm), Standard
505500D	Resolve Open Back Halo • 23-26" (59-66cm), Large
510300D	Resolve Open Back Halo, Ring Only • 17-20" (43-51cm), Small
510400D	Resolve Open Back Halo, Ring Only • 20-23" (51-59cm), Standard
510500D	Resolve Open Back Halo, Ring Only • 23-26" (59-66cm), Large
Halo Vests	
545200C	Halo Vest Assembly • 20-26" (50-66cm), X-Small
545300C	Halo Vest Assembly • 26-33" (66-84cm), Small
545350C	Halo Vest Assembly • 32-36" (81-92cm), Small/Medium
545400C	Halo Vest Assembly • 33-40" (84-102cm), Medium
545450C	Halo Vest Assembly • 38-42" (96-107cm), Medium/Large
545500C	Halo Vest Assembly • 40-46" (102-117cm), Large
545600C	Halo Vest Assembly • 46-52" (117-132cm), X-Large
Lil' Angel Pediatric Vests	
545100K	Lil' Angel Vest, P2 w/ Kodel • P2 (2-6 Years)
545100L	Lil' Angel Vest, P2 w/ Lambswool • P2 (2-6 Years)
545100S	Lil' Angel Vest, P2 w/ Sorbatex • P2 (2-6 Years)
545150K	Lil' Angel Vest, P3 w/ Kodel • P3 (6-12 Years)
545150L	Lil' Angel Vest, P3 w/ Lambswool • P3 (6-12 Years)
545150S	Lil' Angel Vest, P3 w/ Sorbatex • P3 (6-12 Years)

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

SHOES, CAST SHOES

Kelley Cast Boot - Rocker Sole - Open Toe

- Heavy duty canvas upper
- Carpeted inner sole to help reduce wear
- D-ring closure near the toe area for a secure fit
- Upper strap stretches for a firm comfortable fit
- Slip resistant, non-scuff rocker sole
- Hook and loop closure
- Available in 7 sizes
- Sold individually or in cases of 6 (*add 'C' to order number*)



Reorder No.	Size	Length	Unit
0814-8801	X-Small	7 1/4"	Each
0814-8802	Small	9 3/8"	Each
0814-8806	Sm/Med	10 3/8"	Each
0814-8803	Medium	11"	Each
0814-8807	Med/Large	12"	Each
0814-8804	Large	12 5/8"	Each
0814-8805	X-Large	13 1/2"	Each

Comfor™ Tot Cast Boot - Rocker Sole - Open Toe

- Juvenile print canvas upper
- Slip resistant rocker sole
- Hook and loop closure



Reorder No.	Size	Length
0814-7298	Toddler	4 1/2"
0814-7299	Toddler	5 1/2"
0814-7300	Toddler	6 3/8"

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

SHOES, CAST SHOES



Available in 11 sizes!



#882 Cast Shoe - Our most popular style upper, made of strong canvas duck. Straps of canvas duck have hook and loop closures for a more secure fit. The rocker sole aids in the mobility of the patient and the open toe and heel allows for more versatility in sizing.

Order#	Size	sole length
882-98	3XS	4 7/8"
882-INF	2XS	5 1/2"
882-000	Todlr	6 3/8"
882-00	XSM	7 1/4"
882-005	XSM/SML	8"
882-01	SML	9 1/2"
882-015	SML/MED	10 1/2"
882-02	MED	11 1/4"
882-025	MED/LRG	12 1/4"
882-03	LRG	12 3/4"
882-04	XLG	13 1/2"

The Darco SlimLine Cast Boot has taken another step forward in providing cast protection in today's active world. Excellent for use with fiberglass casts or bulky compression bandages because the stretchable upper material conforms for a better fit.

Special Features:



- Square toe design acts as a protective bumper and gives the patient more room to accommodate larger casts and bandages.
- Rocker sole allows for smooth roll through the gait pattern.
- Removable EVA sole provides 10mm of padding to cushion the foot, adding to patient comfort.
- Strapless forefoot closure is easily adjusted to accommodate cast/bandage size and can be switched between left/right closure for ease of re-application.
- Elastic center strap secures neck of the ankle and seats the cast in the shoe.
- Higher ankle strap adds stability

Item No.	DESCRIPTION	SIZE	COLOR	WOMEN SIZE	MEN SIZE
SLQ0B	SlimLine Cast Boot	X-Small	Black	3 - 5	
SLQ1B	SlimLine Cast Boot	Small	Black	5.5 - 8	
SLQ2B	SlimLine Cast Boot	Medium	Black	8.5 - 10	6 - 8
SLQ3B	SlimLine Cast Boot	Large	Black	10.5 - 13	8.5 - 11
SLQ4B	SlimLine Cast Boot	X-Large	Black	13.5 +	11.5 - 14

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

SHOES, POST-OP



Darco's Med-Surg Shoe's strapless forefoot closure eliminates buckle pressure on the forefoot and can expand to accommodate dressings up to 17" in circumference. Available accessories; removable peg insole and toe cover.



ITEM No. (WOMEN)	ITEM No. (MEN)	DESCRIPTION	SIZE	COLOR	WOMEN SIZE	MEN SIZE
MQPB		Med-Surg Shoe	Pediatric	Black		
MQW1B	MQM1B	Med-Surg Shoe	Small	Black	4 - 6	6 - 8
MQW2B	MQM2B	Med-Surg Shoe	Medium	Black	6.5 - 8	8.5 - 10
MQW3B	MQM3B	Med-Surg Shoe	Large	Black	8.5 - 10	10.5 - 12
	MQM4B	Med-Surg Shoe	X-Large	Black		12.5 - 14



ProCare Med/Surg Shoe

Designed for a variety of post-op and post-trauma applications. Foam mesh upper with semi-flexible skid resistant sole and contoured heel. Contact closure straps. Fits right or left foot.

Size	Order # (women)	Womens shoe size	Order # (men)	Mens shoe size
X-Small	79-81142	2 - 4	79-81132	5 - 7
Small	79-81143	4 - 6	79-81133	7 - 9
Medium	79-81145	6 - 8	79-81135	9 - 11
Large	79-81147	8 - 10	79-81137	11 - 13
X-Large	N/A	N/A	79-81138	13+



ProCare Post-Op Shoe

Constructed with E.V.A./rubber, **no-flex** rocker design sole and contoured heel. Low-profile foam/nylon mesh upper has contact closure. Protective tongue provides maximum comfort.

Size	Order # (women)	Womens shoe size	Order # (men)	Mens shoe size
X-Small	79-90192	2 - 4	79-90182	5 - 7
Small	79-90193	4 - 6	79-90183	7 - 9
Medium	79-90195	6 - 8	79-90185	9 - 11
Large	79-90197	8 - 10	79-90187	11 - 13
X-Large	N/A	N/A	79-90188	13+



Baby Post-Op Shoe

Clown print canvas upper, hook and loop closure with a skid resistant sole.

SIZING	CAT. NO.	LENGTH (in)	WIDTH (in)	CHILD'S SHOE SIZE
XX-SMALL	399V	6 1/2	2 1/2	7 1/2 - 10 1/2

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

SHOULDER SUPPORTS

Arco™ Shoulder Support

Dynamic orthosis allows protected shoulder motion.

Provides support for painful, weak, unstable or hypermobile shoulders. Ideal for orthopedic and neurological conditions. Can be worn under clothing. Adjust the three non-stretch straps to restrict shoulder motions such as retraction of the scapula, while limiting rotation and abduction for subluxation conditions. Straps require minimal tightening so the risk of circulatory problems is decreased. Chest strap is padded under the axilla for added comfort. Sleeve opening with "grip" liner allows better humeral fit and less migration. Soft, padded material balances temperature extremes. Hand or machine wash and air dry. Latex free. Fits either the left or right arm.



Sports Shoulder Support by Thermoskin®

The sports shoulder support by Thermoskin provides protection, heat therapy and support for the shoulder. Keeps shoulder warm and flexible while playing sports. Not to be used for dislocations or bracing. Universal design fits either or right or left shoulder.



Added Benefit of Heat Therapy


Only Thermoskin Thermal Supports have been clinically proven to increase subcutaneous skin temperature 2 to 3 degrees for the heat therapy. Thermoskin's patented Trioxon® lining creates a micro-climate that maintains elevated skin temperature while still allowing the skin to ventilate for long term user comfort. The spiral structure of the Trioxon® lining lifts moisture away from the skin and traps air within the lining to help prevent excessive perspiration.

OmoTrain®

The OmoTrain is an active support for treatment of the shoulder. The anatomically-contoured knit ensures optimal adaptation to the shape of the body for maximum comfort. A viscoelastic shaped insert with Velcro tabs can be optionally inserted and positioned as required.

Indications: Shoulder pain (posttraumatic/postoperative, Osteoarthritis of the shoulder)

Mode of action: The support centers the glenohumeral joint and improves function. During movement, the shaped insert produces a local massage effect on the soft tissues of the shoulder.



Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

SPECIALTY SHOES , *Pressure Relief*

The **Darco Peg-Assist System** is built around a multi-purpose removable peg insole that effectively off-loads the plantar aspect of the foot.



Special Features:

- Pegs can be removed to off-load pressure from any area of the plantar surface of the foot.
- Constructed of an 18mm combination of Plastazote®, Poron®, and EVA to provide the perfect combination of shock absorption, off-loading and moldability.
- Peg Support Board adhere to black side of insole to keep pegs from collapsing and keep the off-loaded area intact.

The **OrthoWedge Shoe** has been clinically proven to effectively reduce weight bearing pressure from the forefoot to promote proper healing.

Special Features:

- Square toe design acts as a bumper to provide additional protection when K-wires are present and provides better left/right fit.
- Promotes healing by reducing weight from the forefoot and may be used for any condition from the metatarsal heads distally in which it is desirable to reduce body weight such as ulcers, infections, trauma and surgery.
- 15 degree dorsi-angle shifts weight back to the midfoot and heel.
- Ankle strap eliminates heel slippage and removable forefoot closure eliminates buckle pressure.
- Wide base for extra stability.
- 10mm removable insole fully customizable with over twice the padding.



Also available in HeelWedge version.



The **IPOS Forefoot Relief Orthosis** relieves the forefoot by avoiding any floor contact of the anterior aspect of the foot through a 10° dorsiflexed inclination of the heel and a shortened weight-bearing sole.

- 10° dorsiflexed inclination of the heel.
- Weight-bearing sole ends in the mid foot area.
- Beveled heel for better heel to toe transfer.
- Three closure straps for optimal fit.
- Inner surface completely padded and without seams.
- Non-slip sole for outside use.



Also available in heel relief version.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

SHOULDER PRODUCTS


Saunders Sully Shoulder Support

Traditional shoulder braces have forced you to sacrifice either stability or function. Now you can have both!


On the outside elastic straps attach with hook to the neoprene at any point, in any direction, and with any amount of force you choose. You can functionally stabilize, assist or restrict movement according to the specific needs of each athlete.

On the inside the perforated, breathable, rubber-like neoprene grips the skin so the support moves with the athlete. It's as if the stabilization straps are attached directly to the body.


The Sully comes complete with a highly adaptable, bifurcated neoprene strap, and two additional 4" neoprene straps. Worn by professional and amateur athletes in all sports! Measure biceps in a non-flexed position.



POSTERIOR INSTABILITIES



MUSCLE STRAINS (PECTORIAL ETC.)



ROTATOR CUFF DECELERATION



ANTERIOR INSTABILITIES



INFERIOR INSTABILITIES



MULTI-DIRECTIONAL INSTABILITIES

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

SPINAL HYPEREXTENSION

ACE™ Brace

Three reasons the ACE Brace is simply the most comfortable, cosmetic hyperextension orthosis on the planet!

1. One size fits all backstrap
2. No sternal pressure - no throat impingement
3. Front closure - eases donning by arthritic hands



ACE-ST

Hinged sternal assembly, hinged pubic assembly, universal backstrap and pad, front closure. Ideal for younger, stronger patient who can tolerate direct sternal pressure and is the most economical of our ACE Braces.



ACE-PV

Pectoral pads on sponge rubber disc with velcro® fastener, hinged pubic assembly, universal backstrap and pad, front closure. Ideal for geriatric osteoporotic patient with prominent bony sternum, where orthotist desires to allow patient to shift pectoral pad position for comfort.



ACE-300

A full featured universally sized hyperextension orthosis with one-size-fits-all back strap with innovative Y-tab and swivels on both sternal and pubic plates. Aluminum plates for rapid, precise fit.



ACE-PS

Pectoral pads with adjustable silicone swivels, hinged pubic assembly, universal backstrap and pad, front closure. Ideal for geriatric osteoporotic patient with prominent bony sternum, where orthotist desires pectoral pads to remain exactly as fitted and is the most cosmetic of our ACE Braces.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

SPINAL HYPEREXTENSION

C.A.S.H. Cruciform Anterior Spinal Hyperextension

- Pre-fabricated hyperextension orthosis
- Wide straps
- Fully adjustable height & width
- Available with horizontal or vertical posterior pad
- Cruciform does not make patient contact when worn



Sternal Pad Model

- #CHO100 Fixed sternal and pubic pads
- #CHO300 Hinged sternal and pubic pads
- #CHO400 Hinged sternal pad/fixed pubic pad



Pectoral Extension Model

- #CHO100P Fixed sternal and pubic pads
- #CHO300P Hinged sternal and pubic pads
- #CHO400P Hinged sternal pad/fixed pubic pad

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

SPLINT ROLL

OCL® Splint Roll®-2 Plaster Splinting System



The OCL® Splint Roll®-2 Custom Splinting System is an improved plaster splinting product. Greater strength than other plaster splints. The DuoSet™ plaster gives faster initial set on outer layer. The foam and flannel covering is comfortable for the patient. The open cell foam prevents splint from migrating and improves air exchange for enhanced patient comfort. The flannel reduces moisture retention; ensures faster set times. Non-woven layer allows plaster crystals to lock in the dimpled areas and form additional reinforcement without added weight. Thermal barrier also doubles as a filtering system to prevent plaster from embedding in foam layer. Continuous roll configuration means less waste; cut what you need. Sizing is easy with different widths available. Surged sides provide added patient comfort. Innovative design facilitates custom splinting applications.

Order#	Size
SR215	2"x20'
SR315	3"x20'
SR415	4"x20'
SR515	5"x20'
SR615	6"x20'

Accessories



Disposable buckets make clean-up easier. Economical. Perfect for immersing products; plaster splint roll, plaster bandage and splints, casting tape, synthetic splint roll.

Cutting Paddle protects surface from accidental scratches from knife blade when cutting plaster roll.

Scissors are 9" razor-edge shears to allow for easy cutting of OCL® PolyLITE® Synthetic Splinting System.

OCL® Cutting Knife comes with 3 steel blades stored inside for cutting plaster splinting roll.



Plaster Splint Roll® Dispenser
- 5 compartment



Templates are an aid for preparing custom splints: simply select pattern, trace, prepare splint, cut and apply to patient.



M-PACT® Elastic Bandage High quality elastic retains its shape. Double length bandages available in 4" and 6" widths.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

TEMPERATURE THERAPY



ST600

DeRoyal® Hot/Cold Therapy Unit

- Automatic temperature controlled unit provides hours of hot or cold therapy
- Digital temperature read-out
- Temperature control reduces risk of nerve damage or frostbite
- Measures, adjusts and regulates motor speed for precision temperature control without manual adjustment
- Portable system suitable for any medical setting and continuous treatment at home
- Lightweight, low-profile design is patient friendly
- Automatic shut-off feature lets user know when to add ice
- Hot therapy range (90° - 120°F); cold therapy range (45° - 60°F)
- Hot/cold therapy blankets available in both sterile and non-sterile versions



Knee - Full,
Attached Straps



Knee/Shoulder



Lumbar



Foot/Ankle



Hand/Wrist

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

TEMPERATURE THERAPY

DeRoyal® Hot/Cold Therapy Blankets

Description	Measurement	Qty	Non-Sterile	Sterile
Arm/Leg	10" x 5"	1 Each	ST5010NS	ST5010ST
Lumbar	10½" x 11½"	1 Each	ST5020NS	ST5020ST
Knee/Shoulder	10½" x 10½"	1 Each*	ST5030NS	ST5030ST
Knee/Shoulder - Full	10½" x 10½"	1 Each*	ST5031NS	ST5031ST
Knee - Full, Attached Straps	10½" x 10½"	1 Each*	ST5032NS	ST5032ST
Shoulder - Full, Attached Straps	10½" x 10½"	1 Each*	ST5033NS	ST5033ST
Hand/Wrist - Double	3½" x 5"	2 Each	ST5040NS	ST5040ST
Hand/Wrist	3½" x 5"	1 Each	ST5050NS	ST5050ST
Hip	15½" x 18"	1 Each	ST5060NS	ST5060ST
Cervical	10" x 21½"	1 Each	ST5070NS	ST5070ST
Foot/Ankle	10" x 11½"	1 Each	ST5080NS	ST5080ST
Facial	5½" x 3¾"	1 Each	ST5090NS	ST5090ST

**Also available 10/case*

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

THUMB ABDUCTION

FREEDOM® Thumb Stabilizer™ Soft and safe for the most sensitive hands. Supports arthritic MP joint and reduces pain. Reduces collateral ligament injury pain related to grasp. Strap is positioned through a D-Ring on an angle that permits thumb to be drawn into opposition to support grasp. AliSoft™ Splinting Material and padded Beta Pile™ straps. Wear comfortably all day long. No pressure points, no edge discomfort, and a fully adjustable thumb channel.



Measure width of hand at MP joint from radial aspect of index finger to ulnar aspect of little finger
Suggested code: L3923.

	Small 2 1/4" – 2 3/4"	Medium 2 3/4" – 3 1/4"	Large 3 1/4" – 3 3/4"
Left	5109LS	5109LM	5109LL
Right	5109RS	5109RM	5109RL

\$35.00

Here are some other thumb abduction splints available from Southeast Medical Products...



3PP Thumb Sling
\$32.50



Maramed TGS
\$25.40



ProCare ComfortFORM
Wrist w/Abducted Thumb
\$26.95



ProCare Wrist/Thumb
Wrap
\$9.90 w/o stay
\$17.83 w/stay

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

WALKING BOOTS

MaxTrax walkers are designed to increase stability following trauma or post-operative procedures. Designed with a wider foot bed for comfort and stability, this low profile boot features a rocker bottom which helps promote natural gait and reduce plantar pressures. Cushioned inner and outer sole provides shock absorption to help improve patient comfort during ambulation. Malleable uprights breakdown if needed for custom fit.

MaxTrax™ Walker

MaxTrax™ Air Walker



MaxTrax Tall	MaxTrax Short	Size	Men	Women	MaxTrax Air Short	MaxTrax Air Tall
79-95322	79-95342	Xsmall	NA	Up to 4	79-95422	79-95412
79-95323	79-95343	Small	up to 5	4.5-6	79-95423	79-95413
79-95325	79-95345	Medium	5.5-10	6.5-11	79-95425	79-95415
79-95327	79-95347	Large	10.5-13.5	11.5-14.5	79-95427	79-95417
79-95328	79-95348	Xlarge	14-17	NA	79-95428	79-95418

Southeast Medical Products, Inc
 5524 Commerce Dr
 Orlando FL 32839
 407.851.5848 (ph) • 407.851.9012 (fax)
 800.330.0890 (ph) • 800.851.9012 (fax)

WALKING BOOTS

Wee Walker Boot

The **Bledsoe Wee Walker Pediatric Boot** is the first off-the-shelf low profile pediatric boot that offers superior comfort along with unsurpassed support. Only two sizes are offered thus reducing inventory cost. These sizes will fit left or right, ages 18 months to 5 years old.



Features include:

- Promotes faster healing through strong support with maximum mobility.
- Designed for kids on the go. It's lighter and easier than conventional casts.
- Unique rocker design simulates a near normal gait.
- Keeps kids comfortable with breathable, washable foam wrap.
- Kid-tough performance from durable powder coated aluminum

Children's Shoe Size	Foot Length	Leg Length	Wee Walker Size	Part Number
up to 8	up to 6"	up to 9"	Small - Medium	AL185003
8-12.5	6"-7.5"	9"-12"	Large - XLarge	AL185007

Indications

This device is indicated for post-injury or post-operative application to grade 1, 2, or 3 ankle sprains, and stable or internally fixed fractures of the foot or ankle. Also indicated for treatment of stable or internally fixed fractures or fractures with bridging callous of the middle to distal tibia or fibula.

Contraindications

This device is contraindicated for unstable fractures and proximal tibial or fibular fractures.

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

WALKING BOOTS



...makes all the difference™



Patent Pending

Evenup™ is a new and innovative device that helps to reduce strain of the knee, hip, and back that is often associated with walking in a fracture walker, cast boot or thick soled wound healing shoe.

An **Evenup™** equalizes the induced limb length discrepancy that can cause stress to the rest of the body.

Easily adjustable from 1cm to 2cm.

Keeping patients ambulatory and compliant is what the **Evenup™** is all about!

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839
407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

WRIST ORTHOSES (WH & WHT)

Wrist Hand Orthosis

Provides proper support and immobilization of the wrist. Allows unobstructed motion of the thumb and fingers.

Suggested L-Code: L3984

Indications

- Post cast fracture immobilization
- Wrist fractures and sprains
- Cumulative Trauma Disorder
- Carpal Tunnel Syndrome
- Post tendon transfer support
- Postoperative support after removal of internal or external fixation devices



Wrist Hand Thumb Orthosis

Provides proper support and immobilization of the wrist and thumb. Allows unobstructed motion of the fingers.

Suggested L-Code: L3984

Indications

- Post cast fracture immobilization
- Wrist fractures and sprains
- Cumulative Trauma Disorder
- Carpal Tunnel Syndrome
- Post tendon transfer support
- Postoperative support after removal of internal or external fixation devices
- Gamekeeper's thumb
- DeQuervain's Syndrome



Size	M-P Diameter*	Wrist Circum.	Forearm Circum.	Length of Brace	Wrist Hand product number	Wrist Hand Thumb product number
Small	2¾ - 3¼	6 - 7	8 - 10	7½	WH100 (L or R)	WHT100 (L or R)
Medium	3¼ - 3¾	7 - 8	10 - 12	7½	WH200 (L or R)	WHT200 (L or R)
Large	3¾ - 4	8 - 9	12 - 14	8	WH300 (L or R)	WHT300 (L or R)

Southeast Medical Products, Inc
 5524 Commerce Dr
 Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
 800.330.0890 (ph) • 800.851.9012 (fax)

WRIST SUPPORTS

ComfortFORM™ Wrist

Durable lightweight foam laminate is Lycra® lined for breathability and patient comfort. Preformed aluminum stay and loop/lock closure helps provide anatomically correct fit and proper support. Ideal for sprains, strains and control of carpal tunnel syndrome symptoms.



Left	Right	Wrist Circ.	Length	Size
79-87291	79-87281	Up to 4.5"	7.0"	Ped/XXS
79-87292	79-87282	4.5"-5.5"	7.5"	X-Small
79-87293	79-87283	5.5"-6.5"	7.5"	Small
79-87295	79-87285	6.5"-7.5"	8"	Medium
79-87297	79-87287	7.5"-8.5"	8.5"	Large
79-87298	79-87288	8.5" +	8.5"	X-Large

ComfortFORM™ Wrist w/Abducted Thumb

Durable lightweight fabric is Lycra® lined for breathability and patient comfort. Dual aluminum stays support palmer surface of wrist and extensor surface of thumb. Stays can be adjusted for proper angulation. Contact closure straps. Ideal for wrist sprains, strains, scaphoid injuries, Carpal Tunnel, Gamekeeper's thumb, and DeQuervain's Syndrome symptoms.



Left	Right	Wrist Circ.	Length	Size
79-87312	79-87302	up to 5.5"	7.5"	X-Small
79-87313	79-87303	5.5"-6.5"	7.5"	Small
79-87315	79-87305	6.5"-7.5"	8"	Medium
79-87317	79-87307	7.5"-8.5"	8.5"	Large
79-87318	79-87308	8.5"+	8.5"	X-Large

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

WRIST SUPPORTS

Premier Wrist Brace

- Indicated when wrist immobilization and full hand function are desired
- Maximum contouring at the palmar crease allows for full finger function
- Additional contouring at the fifth M.P. allows for palmar spread
- Removable and adjustable palmar stay for maximum support
- Generous thumb opening permits freedom of movement
- Constructed of durable, breathable, padded material for maximum wear and comfort
- Loop and lock closure for easy adjustment
- 8" length



Reorder No.	Reorder No.	Wrist Circumference
Left	Right	
0814 4390	0814 4380	XX-Small (up to 5 3/4")
0814 4391	0814 4381	X-Small (5 3/4" - 6 1/2")
0814 4392	0814 4382	Small (6 1/2" - 7 1/4")
0814 4393	0814 4383	Medium (7 1/4" - 8 1/4")
0814 4394	0814 4384	Large (8 1/4" - 9")
0814 4395	0814 4385	X-Large (9" - 10")

Premier Wrist Brace with Thumb Spica

- Indicated when immobilization of the MP and CMC joints is desired
- Indicated for deQuervain's syndrome and rheumatoid arthritis
- Maximum contouring at the palmar crease allows for full finger function
- Additional contouring at the fifth M.P. allows for palmar spread
- Removable and adjustable palmar stay for maximum support
- Removable rigid thumb stay is adjustable for a custom fit
- Constructed of durable, breathable, padded material for maximum wear and comfort
- Loop and lock closure for easy adjustment
- 8" length



Reorder No.	Reorder No.	Wrist Circumference
Left	Right	
0814 4551	0814 4561	X-Small (up to 6 1/2")
0814 4552	0814 4562	Small (6 1/2" - 7 1/4")
0814 4553	0814 4563	Medium (7 1/4" - 8 1/4")
0814 4554	0814 4564	Large (8 1/4" - 9")
0814 4555	0814 4565	X-Large (9" - 10")

Southeast Medical Products, Inc
5524 Commerce Dr
Orlando FL 32839

407.851.5848 (ph) • 407.851.9012 (fax)
800.330.0890 (ph) • 800.851.9012 (fax)

How to Order

By Phone (M-F, 9am-5pm) 407.851.5848 or 800.330.0890
By Fax (24/7/365) 407.851.9012 or 800.851.9012
By Email (24/7/365) orders@southeastmedical.com

Payment Terms

Payment terms are Net 30 days. For your convenience, we accept payment by Visa, MasterCard and American Express.

Warranty Policy

Products sold by Southeast Medical Products are covered by the manufacturer's warranty, which varies by manufacturer. Most products are generally warranted against defects in materials and workmanship for a period of at least six months from the date of purchase.

Return Policy

All returns must be authorized by Southeast Medical Products. Before returning any product for any reason, please call 1.800.330.0890 to obtain a Return Authorization Number. Non-warranty returns will be accepted within 30 days from date of purchase and may be subject to a restocking fee.